




These Safe Basics general purpose sponges are highly absorbent and offer an excellent alternative to cotton gauze sponges.

- Soft, 30 gram weight viscose/poly blend, non-woven material for increased patient comfort.
- No chemical additives to avoid allergic reactions.
- Latex free. Non-sterile.
- Can be used for patients, disinfecting and general cleaning.







Use fewer sponges per extraction with highly absorbent, cotton-filled gauze.

- 100% white cotton with cotton-fibre filling.
- Latex free, multi-purpose sponges save time.
- Non-sterile.

Item No.	Description	Unit/Pkg	Pkg/Case
4521	2x2 N/S Dental Cotton Gauze 8 Ply	200	25





A wide variety of sizes and ply options allows you to select the appropriate option for different procedures.

- All-gauze sponges with high absorbency allow for the use of fewer sponges resulting in savings of time and money.
- 100% pure, white cotton.
- Non-sterile.

Description	Unit/Pkg	Pkg/Case
2x2 N/S Cotton Gauze 8 Ply	200	25
2x2 N/S Cotton Gauze 12 Ply	200	40
3x3 N/S Cotton Gauze 8 Ply	200	20
3x3 N/S Cotton Gauze 12 Ply	200	20
4x4 N/S Cotton Gauze 8 Ply	200	20
4x4 N/S Cotton Gauze 12 Ply	200	10
8x4 N/S Cotton Gauze 12 Ply	200	10
	2x2 N/S Cotton Gauze 8 Ply 2x2 N/S Cotton Gauze 12 Ply 3x3 N/S Cotton Gauze 8 Ply 3x3 N/S Cotton Gauze 12 Ply 4x4 N/S Cotton Gauze 8 Ply 4x4 N/S Cotton Gauze 12 Ply	2x2 N/S Cotton Gauze 8 Ply       200         2x2 N/S Cotton Gauze 12 Ply       200         3x3 N/S Cotton Gauze 8 Ply       200         3x3 N/S Cotton Gauze 12 Ply       200         4x4 N/S Cotton Gauze 8 Ply       200         4x4 N/S Cotton Gauze 12 Ply       200





Peace of mind with sterile, heat-sealed and secure envelope for added safety from cross contamination where additional caution is warranted.

- Absorbent, soft, non-linting and breathable material for increased patient comfort.
- Convenient, per-use, individual peel-open design.
- 2 sponges per envelope.

Item No.	Description	Unit/Pkg	Pkg/Case
2240	2x2 N/W Sterile Sponge 4 Ply	50	30
3340	3x3 N/W Sterile Sponge 4 Ply	50	24
4440	4x4 N/W Sterile Sponge 4 Ply	50	12





Increased convenience with easy-to-use, individual, peelopen package design.

- Highly absorbent, soft sponges for optimal protection and patient comfort.
- Sterile, heat-sealed envelope for superior safety.
- 100% cotton.
- 1 sponge per envelope.





Increased convenience with easy-to-use, individual, peel-open package design.

- Highly absorbent, soft sponges for optimal protection and patient comfort.
- Sterile, heat-sealed envelope for superior safety.
- 2 sponges per envelope.

Item No.	Description	Unit/Pkg	Pkg/Case
2212	2x2 Woven Sterile Sponge Single 12 Ply	100	24
3312	3x3 Woven Sterile Sponge Single 12 Ply	100	24
4412	4x4 Woven Sterile Sponge Single 12 Ply	100	12

Item No.	Description	Unit/Pkg	Pkg/Case
2120	2x2 Woven Sterile Sponge Double 12 Ply	100	30
2800	2x2 Woven Sterile Sponge Double 8 Ply	100	30
3120	3x3 Woven Sterile Sponge 12 Ply	80	30
3800	3x3 Woven Sterile Sponge Double 8 Ply	80	30
4120	4x4 Woven Sterile Sponge Double 12 Ply	50	24
4800	4x4 Woven Sterile Sponge Double 8 Ply	50	12