

**Distech Cotton Dental Rolls Non Sterile Medium #2**

---

|                             |                                                          |
|-----------------------------|----------------------------------------------------------|
| <b>Item Code</b>            | 3554                                                     |
| <b>Country of Origin</b>    | China                                                    |
| <b>Primary Material</b>     | 100% Cotton                                              |
| <b>Length</b>               | 38 mm                                                    |
| <b>Diameter</b>             | 10 mm                                                    |
| <b>Retention (XoW)</b>      | > 4                                                      |
| <b>Absorption (seconds)</b> | < 30                                                     |
| <b>Sterility</b>            | Non-Sterile                                              |
| <b>Color</b>                | White                                                    |
| <b>Usage</b>                | For Single-use only                                      |
| <b>Packaging</b>            | 2000 rolls / bag<br>12 bags / case = 24,000 rolls / case |



**This product is 100% Latex-Free**