

TECHNICAL DATA SHEET

Cotton-Filled Exodontia Sponge Item # 4521

Intended Use

An absorbent sponge for external use to be placed directly on a patient's wound to absorb exudates.

Characteristic	Specifications		
Description	Folded gauze filled with cotton		
Material composition	100% Cotton. Latex free		
Color	White		
Ply	8 Ply		
Dimensions	Size	Width (cm)	Length (cm)
	2" x 2"	5,0	5,0
Packaging	Sponges/Bag	Bags/Case	Sponges/Case
	200	25	5000

Item Code	Description	Color	Specific Attributes
4521	Medicom Cotton-Filled Exodontia Sponges, 8 ply, 2" x 2", 200 Sponges/Bag, 25 Bags/Case, 5000 Sponges/Case	White	Non-Sterile

Regulatory

Health Canada Class 1
FDA Class 1, Product Code NAB

Visual

Inner Sleeve

Storage Conditions

Store under cool dry place.

Country of Origin

China

