

### TECHNICAL DATA SHEET

# SafeBasics Non-Woven Sponge Items # 2100-CH to 2103-CH

#### **Intended Use**

A hydrophilic non-sterile sponge intended to cover a wound, such as to be placed directly on a patient's wound to absorb exudates.

Characteristic	Specification				
Description	Soft non-woven material, folded				
Material Composition	30% Polyester / 70% Viscose				
Ply	4 ply				
Color	White				
Dimension	Size	Item Code	Width (cm)	Length (cm)	
	2" x2"	2101-CH	5,0	5,0	
	3" x 3"	2103-CH	7,5	7,5	
	4" x 4"	2100-CH	10,0	10,0	
Packaging	Item Code	Pouches/ Sleeve	Sleeves/ Case	Sponges/ Case	
	2101-CH	200	25	5000	
	2103-CH		20	4000	
	2100-CH		10	2000	

Item Code	Description	Color	Specific Attribute
2100-CH	Medicom SafeBasics General Purpose Non-Woven Sponges, 4 ply, 4" x 4" (10 x 10 cm), 200 Sponges/Sleeve, 10 Sleeves/Case, 2000 Sponges/Case	White	Non-Sterile
2101-CH	Medicom SafeBasics General Purpose Non-Woven Sponges, 4 ply, 2" x 2" (5 x 5 cm), 200 Sponges/Sleeve, 25 Sleeves/Case, 5000 Sponges/Case	White	Non-Sterile
2103-CH	Medicom SafeBasics General Purpose Non-Woven Sponges, 4 ply, 3" x 3" (7,5 x 7,5 cm), 200 Sponges/Sleeve, 20 Sleeves/Case, 4000 Sponges/Case	White	Non-Sterile

#### Visual

Inner Sleeves



## Regulatory

Health Canada Class 1 FDA Class 1, Product Code NAC

### **Storage Conditions**

Store in a dry, cool place

# **Country of Origin**

China

Marketing: TDS 3.4.4-3 Date: 27-Nov-2015 Page 1 of 1