



TECHNICAL DATA SHEET
Non-Woven Sponges
Items #A2336 to A2338 & A7632 to A7634

Intended Use

May be used for one or more of the following purposes: cleansing wounds or skin, absorbing biological fluids or support organs and tissue during surgical procedures

Characteristics	Specification			
Description	Soft non-woven sponges			
Material	Viscose: 70% Polyester: 30% Latex-free			
Color	White			
Ply	4 ply			
Dimension ⁽³⁾	Item Code	Width (cm)	Length (cm)	
	A2336	20.4 (8")	10.2 (4")	
	A2337	10.2 (4")	10.2 (4")	
	A2338	5.1 (2")	5.1 (2")	
	A2339	7.6 (3")	7.6 (3")	
	A7632	5.1 (2")	5.1 (2")	
	A7633	7.6 (3")	7.6 (3")	
	A7634	10.2 (4")	10.2 (4")	
Sterilization	Ethylene Oxide			
Packaging	Item Code	Qty/Pouch	Tray or Bag/Case	Qty/Case
	A2336	2	20	1000
	A2337		24	1200
	A2338		60	3000
	A2339		48	2400
	A7632	200/Bag	40	8000
	A7633		20	4000
	A7634		10	2000

Storage Conditions

Store under room temperature and dry conditions