

Quala Dental Products

small equipment



5200 Sweep with Timer Ultrasonic Cleaning Unit

Quala Ultrasonic Cleaners are manufactured with vinyl-clad housing, which ease cleaning and resist scratches. All Quala machines include a stainless steel tank and cover, reducing the risk of cross-contamination and are equipped with a 60-minute timer. All Quala Ultrasonic Cleaners are backed by a 2½ year warranty.

Quala Sweep Machines utilize sweep technology, a cleaning wave that sweeps back and forth throughout the ultrasonic tank, reaching all corners of the machine – no reduced cleaning areas.

FEATURES

Sweep cleaning

Swivel drain

Cover included

BENEFITS

Sweep cleaning wave reaches all corners of the machine; thorough and complete cleaning without weak cleaning areas

Convenient system of draining the machine

Reduces the possibility of water splashing to work area

CODE	DESCRIPTION
5200 Sweep with Timer	
1209	Internal tank dimensions: 11¾" x 6" x 4" (29.8 x 15.2 x 10.2 cm) Capacity: 4¼ quarts (4 liters)
5200 Accessory Kit	
1212	Includes one solid side basket, positioning cover, two 600 mL beakers complete, one bur tray and one bur box

SWEEPZONE Performance Power

Automatically controlled frequency changes create an ultrasonic cleaning "grid" which literally "sweeps" through the entire tank, unlike that pattern created by traditional ultrasonic systems. The result is total cleaning, which means constant power and consistent output.



Traditional Ultrasonic System



SweepZone Ultrasonic System

Diminished cleaning power
Lowest cleaning power

Indicates areas of constant, intense cleaning power

Distributed by: