



HANDLER

RED WING INTERNATIONAL



Your **New Lab** Design Starts Here!

We offer full lab functionality and performance to fit in any space and any budget

LAB FURNITURE



AMERICAN-MADE DURABILITY IN DENTAL LAB EQUIPMENT

HANDLER MANUFACTURING, INC.
www.handlermfg.com

Toll Free: 800.274.2635

612 North Avenue E., Westfield N.J. 07090
Tel (908) 233-7796 Fax (908) 233-7340

Our Color Selections:



Warm Gray, Colonial White,
Designer Beige, Eagle White,
Blue Heaven, Charcoal, Merlot



Modular, Flexible Design and Easy Assembly

Custom Benches and Cabinets
Built to Your Specifications.
Call us for a free consultation.




All Benches and Cabinets Available in:


2-TECHNICIAN CORNER WORK STATION



311 BENCH



253 CORNER CABINET



213 CERAMIST BENCH

Benches Specially Designed for Dental Lab Operations


SMALL FULL SERVICE LAB






DENTIST OFFICE LAB

Modular Cabinets for Flexible Design & Ease of Assembly





GROUP PRACTICE




PTL-215SP BENCH


Customized Solutions

- Custom Benches and Cabinets Built to Your Specifications
- Stainless Steel Tops with Raised Edge to Prevent Water Runoff
- Special Shelving for Storage or Holding Burnout Ovens
- Slide Out Trays Installed for Cabinet Curing Units or Other Needs

4-TECHNICIAN DENTURE WORK BENCH



215WS BENCH



MULTI-STATION TECH RITE



211 BENCH

Durable, Rugged
Resistop[®] Countertops


All Benches and Cabinets Available in:




CORNER WORK STATION



Model 213CWS



231 BENCH



213 BENCH



COMMERCIAL - HOSPITALS, GOVERNMENT AND EDUCATION



PRO TECH *Lite*

Sensible Design Combining Art & Engineering

Now full lab functionality and performance fits comfortably in any space and any budget.

The PTL-215SP is engineered to provide affordable and durable functionality with style and performance. As the newest addition to Handler's Pro Tech Lite bench line, the PTL-215SP is the most comprehensive work station there is. Standard with the PTL-215SP is a full 24" deep work surface free of obstacles for maximum work area and the list of standard components is long.

- All 18 gauge welded steel construction with an oven-cured powder-coated painted finish
- Spacious 6" riser or chase to run air, gas and electric lines
- Totally enclosed, powerful and quiet 62-II Dust Collector pre-wired for ease of operation
- Double tier shelf to hold cases conveniently
- Dual fluorescent bulb overhead lighting mounted on stylish and strong support with the pre-wired light switch
- Six outlet pre wired power strip
- Removable 16" x 12" stainless steel work surface protector
- Laboratory grade solid, smooth Resistop® work surface that resists stain and scratches
- Two suction inlet options
 - Bench top 3" suction inlet with blast gate
 - In drawer suction complete with protective acrylic shield and padded comfort arm rest
- Four drawers for storing tools and equipment
- Easy set-up and assembly
- Dimensions: 42"W x 30"D x 32" H work surface



PTL-215SP



Streamline your lab operation with the **New Handler Pro Tech Lite Bench Series**

PTL-211 Workstation

Shipping weight: 175lbs
Dimensions: 36" W x 34" H x 26" D



- Sturdy steel construction
- Leg levelers
- Durable oven-cured powder-coated finish
- Light weight – 34"W x 24"D Work surface
- Power strip
- In drawer suction
- Resistop© work surface with surround wrap
- Stainless steel work surface protector
- Pre-installed air hose and rubber blow gun
- Padded, removable comfort arm rest
- Compact space design
- European styling
- No assembly – easy install
- Affordable price

PTL-231 Petite Plaster Bench

Shipping weight: 185lbs
Dimensions: 35" W x 36" H x 24" D



- Sturdy steel construction
- Leg levelers
- Durable oven-cured powder coated finish
- Light weight – 35"W x 24"D work surface
- Resistop© sink and drain
- Gooseneck lab faucet
- Optional in drawer plaster bin liner
- Optional kick plate hands free water activation
- Resistop© work surface
- Compact space design
- No assembly – easy install
- Affordable price

www.handlermfg.com
800.274.2635

Contact Handler for Free lab layout service.
No lab is too small or too big.

LABORATORY WORK BENCHES

Handler's work stations and cabinets are built of heavy gauge sheet steel and powder coated for longevity and durability.

211 TECH-RITE BENCH

Five drawer technician's bench with resistop. Resistop is similar to man made slate. Built of heavy gauge sheet steel with heavy gauge tubular steel legs, the 211 is shipped complete with four legs with leg levelers, duplex electrical outlet, suction inlet and knock-outs for air and gas valves. This bench has an oven-cured, powder-coated finish and Resistop® Countertops. *NOTE: Gas and air valves as shown are at additional charge.

Shipping Wgt.: 170 lbs. [77.1 Kg]

Unit Dimensions: 42" W x 24 1/2" D x 32" to work surface
[107cm x 62cm x 81cm]

with 42" W x 6 1/2" D x 6" H

Riser [107cm x 17cm x 15cm], mounted rear



213 CERAMIST BENCH

This seven drawer bench is designed particularly for the ceramist. In addition to the two larger side drawers which will accommodate large porcelain jars, its side cabinet contains a hinged door and compartment for a vacuum pump or other accessories. The side cabinet is shipped for assembly on the right side of the bench. Built of heavy gauge sheet steel, the Tech-Rite Ceramist Bench is shipped complete with legs, duplex electrical receptacle, suction inlet and knock-outs for air and gas valves. This bench has an oven-cured, powder-coated finish and Resistop® Countertops.

Specify placement of side cabinet - LEFT or RIGHT.

Shipping Wgt.: 270 lbs. [122.4 Kg]

Unit Dimensions: 58" W x 24 1/2" D x 32" to work surface
[147cm x 62cm x 81cm].

Bench contains 58" W x 6 1/2" D x 6" H

Riser [147cm x 17cm x 15.2] mounted rear

215SP SINGLE PEDESTAL BENCH

Five drawer single technician's bench with Resistop tops. Bench is supported by a 17" deep x 42" wide [28cm x 107cm] single pedestal. The 215SP contains a 42" wide x 6 1/2" deep x 6" high riser [107cm x 17cm x 15cm] mounted rear with one duplex receptacle, suction inlet and knock outs for air and gas valves in the vertical riser. This bench must be secured to the wall of the laboratory or installed to the floor as a free standing bench. May be combined with additional benches to form a free standing or wall mounted island of benches. Constructed of heavy gauge sheet steel with an oven-cured, powder-coated finish and Resistop® Countertops.

Shipping Wgt.: 225 lbs. [102 Kg]

Unit Dimensions: 42" W x 30 1/2" D x 32" to work surface
[107cm x 77cm x 81cm]





215WS PEDESTAL BENCH

A two technician pedestal bench containing five drawers for each technician. Benches are mounted face to face supported by pedestal. Bench also includes a 60U dust collector for two technicians, two fishmouths, flexible ducting, clamps, blastgates and two lighted remote switches mounted on vertical risers.

Supplied with Resistop® tops, and an oven-cured, powder-coated finish. *NOTE: Gas and air valves as shown are at additional charge.

Shipping Wgt.: 450 lbs. [204.1 Kg]

Unit Dimensions: 42" W x 49" D x 32" to work surface [107cm x 125cm x 81cm]

Riser is 42"W x 13" D x 6" H [107cm x 33cm x 15cm]

Available without dust collector - Model 215

800 DEXTERITY 2000 WORKSTATION

Handler has developed a new workstation to satisfy the ergonomic requirements of the laboratory technicians. The Dexterity 2000 Workstation is a modern, practically designed bench which is sturdy, comfortable, and durable. Workstations available in 42", 48" and 60" wide. Varying widths are available on special order.

800

Mechanically Adjustable Height Bench 42"W X 30"D Steel construction with laminated top workstation. 12"D X 30"W shelf, dust collector, light, two duplex outlets and knock outs for gas and air. Shipping Wgt: 500 lbs

801

Mechanically adjustable height bench 48"W x 30"D with laminated workstation top. 12"D X 30"W shelf, dust collector, light, two duplex outlets and knock outs for gas and air. Shipping Wgt: 525 lbs

802

Mechanically adjustable height bench 60"W x 30"D with laminated workstation top. 12"D X 30"W shelf, dust collector, light, two duplex outlets and knock outs for gas and air. Shipping Wgt: 575 lbs

P901-04: EXTRA SHELF, 12"D X 30" W SOLD SEPARATELY

P800-1: SET OF THREE DRAWERS SOLD SEPARATELY



*Shown as Model 802 with extra shelf (P901-04) and drawers (P800-1)



311 MASTER TECHNICIAN BENCH

The ultimate bench for the accomplished technician, designed to provide maximum work space for supervision and teaching. This bench contains 13 drawers, a 36" [91cm] electrical strip, air and gas valves, and a 3" [7.6cm] suction inlet with blastgate. This bench is constructed of heavy gauge steel with an oven cured, powder coated finish and includes a Resistop work surface. A 6 1/2" deep x 6" high [17cm x 15cm] vertical riser is mounted in the rear of the bench.

Shipping Wgt.: 350 lbs. [158.7Kg]

Unit Dimensions: 74" W x 24 1/2" D x 32" H to work surface [188cm x 62cm x 81cm]



231 PLASTER BENCH

A seven foot bench built for use as a plaster bench. Bench contains one 35" [89cm] double door cabinet, one 24" [61cm] single door cabinet and one 24" [61cm] four drawer cabinet combined. Built-in 16" x 12" x 8" deep Resistop® sink with chrome plated mixing faucet, and adjacent plaster (waste) chute. Bench constructed of heavy gauge sheet steel with an oven-cured, powder-coated finish with Resistop® Countertop and backsplash.

Shipping Wgt.: 475 lbs. [215.4 kg]

Dimensions: 83 3/8" W x 24 1/2" D x 36 3/4" H to work surface
[213cm x 62cm x 109cm]

Available with 6" riser at additional cost.

247DL – In-drawer liner available for plaster and stone storage.

231HF – Optional hands-free water faucet activator.



235 CASTING/SOLDERING BENCH

A seven foot long bench designed for casting and soldering. This bench consists of one 35" [90cm] double door cabinet, one 24" [61cm] single door cabinet and one 24" [61cm] single door cabinet with tool area insulated for placement of torches and tools. This cabinet's work surface is covered with Resistop material and contains a cut-out with a 19" Diameter x 14 1/2" High casting well.

Available at additional cost with 6" [15cm] high riser in lieu of backsplash making bench 29" [73cm] in depth. Constructed of heavy gauge sheet steel with an oven-cured, powder-coated finish.

Shipping Wgt.: 460 lbs. [208.6 kg]

Dimensions: 83 3/8" W x 24 1/2" D x 36 3/4" H to work surface
[213cm x 62cm x 109cm]



235CS CASTING/SOLDERING BENCH

All the benefits of the seven foot long 235 Casting & Soldering bench in a compact, ergonomic design. The unit comes complete with a cutout and casting well, 3 drawers for storage, knock-outs for air & gas connection, a 6" riser and two duplex outlets. Available with a backsplash in place of the 6" riser; bench depth will be 23".

Shipping Weight: 340 lbs. [155 kg]

Dimensions: 48" W x 29" D x 36 3/4" H [122cm x 74cm x 109cm]

LABORATORY CABINETS & BENCHES

241 DOUBLE DOOR CABINET

A 35" [90cm] wide double door cabinet built of heavy gauge sheet steel and covered with Resistop work surface. Full size doors provide complete access to enclosed area. Cabinets and doors finished in an oven-cured, powder-coated finish.

Shipped with 6" Backsplash, may be ordered with 6" Riser option at additional cost.

Shipping Wgt.: 185 lbs. [83.9 kg]

Dimensions: 35" W x 23 3/4" D x 36 3/4" H
[89cm x 58cm 92cm]



243 SINGLE DOOR CABINET

A 24" [61cm] wide single door cabinet built of heavy gauge sheet steel with Resistop top. Full sized door allows complete access to inside of cabinet.

Shipped with 6" Backsplash, may be ordered with 6" Riser option at additional cost.

Shipping Wgt.: 130 lbs. [59 kg]

Dimensions: 24 1/8" W x 23 3/4" D x 36 3/4" H
[61cm x 63cm x 92cm]

245 SINGLE DOOR CABINET

An 18" [46cm] wide single door cabinet constructed of heavy gauge sheet steel with Resistop work surface. Complete access to cabinet interior is afforded by full size door. Shipped with 6" Backsplash, may be ordered with 6" Riser option at additional cost. Cabinet finished in an oven-cured, powder-coated finish.

Shipping Wgt.: 115 lbs. [52.2 kg]

Dimensions: 18 1/8" W x 23 3/4" D x 36 3/4" H
[46cm x 76cm x 92cm]



247 FOUR DRAWER CABINET

This four drawer cabinet is manufactured of heavy gauge sheet steel and finished in a durable, baked coating. Each drawer is equal in size. Work surface covered in Resistop material. Shipped with 6" Backsplash, may be ordered with 6" Riser option at additional cost.

Shipping Wgt.: 160 lbs. [72.6 kg]

Dimensions: 24 1/8" in W x 23 3/4" D x 36 3/4" H [61cm x 58cm x 92cm]



251 CORNER CABINET

A single door cabinet designed to be placed in a corner and abut modular benches on either side. Constructed of heavy gauge sheet steel and finished in an oven-cured, powder-coated finish. Shipped with 6" Backsplash, may be ordered with 6" Riser option at additional cost.

Shipping Wgt.: 190 lbs. [86.2 kg]

Dimensions:

Wall side - 46" [89cm], abutment side - 23 1/2" [60cm],

Door side - 18 1/8" [46cm],

Height - 36 3/4" [93cm]



247DL - In-drawer liner/separator available to hold plaster and stone.



256 PACKING BENCH

A packing bench utilized for general lab procedures. Manufactured of heavy gauge sheet steel with an oven-cured, powder-coated finish, this bench has a myriad of uses in the laboratory.

It is supplied with a Resistop work surface, a 6" vertical riser and two duplex electrical receptacles.

Shipping Wgt.: 145 lbs. [65.8 kg]

Dimensions: 42" W x 24 1/2" D x 32" H to work surface
[107cm x 62cm x 81cm]

261 BOILOUT CABINET

Boilout bench is constructed of heavy gauge sheet steel with added structural supports to insure the secure containment of your boilout unit. Bench contains two front doors with louvers to allow for proper air circulation throughout heated boilout unit. The piece is coated in a durable, baked finish. Shipped with 6" Backsplash, may be ordered with 6" Riser option at additional cost.

Shipping Wgt.: 215 lbs [97.5 Kg]


Dimensions: 35" W x 23 3/4" D x 36 3/4" H [89cm x 58cm x 92cm]

NOTE: Shown with 26102 Boilout/Curing unit—NOT INCLUDED. Specify cut out size for Boil-Out Unit.

Contact Handler prior to ordering.



**Shown with 26102
Dual Tank Curing System**



275 WALL CABINET

A double sliding door cabinet built to be mounted on the wall. Cabinet manufactured of heavy gauge sheet steel with adjustable single interior shelf. Finished in an oven-cured, powder-coated finish with glide-rite doors.

NOTE: Additional shelves available at additional cost.

Shipping Wgt.: 100 lbs [45.4 Kg]

Dimensions: 36" W x 10" D x 30" H [91cm x 25cm x 76cm]



222 DOUBLE SIDE WORK SHELF

Designed for placement on top of riser where benches are placed face to face. The shelf is constructed of heavy sheet steel and perforated steel. Built to hold case pans and finished in a baked, durable coating. Specify color when ordering.

Shipping Wgt.: 75 lbs [34 Kg]
Dimensions: 42" W x 15" D x 12" H
[107cm x 38cm x 30cm]



223 SINGLE SIDE WORK SHELF

Designed for placement on top of single position work benches. Constructed of heavy sheet steel and perforated steel, the shelf attaches to riser. Finished in a durable coating. Specify color when ordering.


Shipping Wgt.: 55 lbs. [24.9 kg]
Dimensions: 42" W x 7 1/2" D x 12" H
[107cm x 19cm x 30cm]



213SDT DOUBLE TIER SINGLE SIDE

Double tier work shelf for ceramist bench, constructed of heavy sheet steel and perforated steel, attaches to riser. Built to hold case pans. Finished in a durable baked coating. NOTE: Contact Handler for pricing on special size work shelves of different lengths and widths. Specify color when ordering.

Shipping Wgt.: 80 lbs. [36.3 kg]
Dimensions: 58" W x 7 1/2" S x 24 " H
[147cm x 19cm x 61cm]



222DT DOUBLE SIDE/DOUBLE TIER WORK SHELF

Designed for placement on top of riser where benches are placed face to face. The shelf is constructed of heavy sheet steel and perforated steel. Built to hold case pans and finished in a baked, durable coating. Second shelf doubles work area. Specify color when ordering.


Shipping Wgt.: 95 lbs [34 Kg]
Dimensions: 42" W x 15" D x 24" H
[107cm x 38cm x 61cm]



223DT DOUBLE TIER SINGLE SIDE

Designed for placement on top of single position work benches. Constructed of heavy sheet steel and perforated steel, the shelf attaches to riser. Finished in a durable coating. Second tier doubles workshelf area. Specify color when ordering.


Shipping Wgt.: 75 lbs. [34 kg]
Dimensions: 42" W x 7 1/2" D x 24" H
[107cm x 19cm x 61cm]



213DT DOUBLE TIER DOUBLE SIDE

Designed for placement on top of riser where benches are placed face to face. The shelf is constructed of heavy sheet steel and perforated steel. Built to hold case pans and finished in a baked, durable coating. Specify color when ordering.

Shipping Wgt.: 90 lbs
Dimensions: 58" W x 15" D x 24" H
[147cm x 38cm x 61cm]



Handler's benches and cabinetry can withstand the rigors of the lab environment and are manufactured of 18 gauge steel with long-lasting Resistop® work surface.

Handler's solutions are the most affordable solutions in the industry. The staff at Handler has the experience and know how for efficient lab work flow and can customize almost any solutions to meet requirements.

No project is too small or too big



Handler manufactures its products in the USA and offers turn around times that are very reasonable.

Handler offers a variety of colors and styles to meet your needs.

Handler provides the necessary accessories such as gas and air valves, boil out and curing units, casting areas, plaster work areas, polishing units, central or decentralized dust collection systems, hands free faucets and sinks, comfortable seating, shelves, lighting, plaster bins, wall mounted storage and more.

The staff at Handler understands labs and we care about our customers!

AMERICAN-MADE DURABILITY IN DENTAL LAB EQUIPMENT

FAMILY OWNED & OPERATED SINCE 1920

- All 18 gauge welded steel construction with an oven-cured powder-coated painted finish
- Spacious 6" riser or chase to run air, gas and electric lines
- Totally enclosed, powerful and quiet 62-II Dust Collector pre-wired for ease of operation
- Double tier shelf to hold cases conveniently
- Dual fluorescent bulb overhead lighting mounted on stylish and strong support with the pre-wired light switch
- Six outlet pre-wired power strip
- Removable 16" x 12" stainless steel work surface protector
- Laboratory grade solid, smooth Resistop® work surface that resists stain and scratches
- Two suction inlet options
 - Bench top 3" suction inlet with blast gate
 - In drawer suction complete with protective acrylic shield and padded comfort arm rest
- Four drawers for storing tools and equipment
- Easy set-up and assembly
- Dimensions: 42"W x 30"D x 32" H work surface



PRO TECH Life


Made in the U.S.A.
Since 1920

HANDLER

RED WING INTERNATIONAL

(800) 274 2635

612 North Ave East, Westfield, NJ 07090

Tel (908) 233-7796 Fax (908) 233-7340

WWW.HANDLERMFG.COM

PTL-215SP