



HANDLER

16B CHUK/CHANGER W/LATHE

The Handler Chuk/Changer was designed to provide full automatic, precision changing of 1/4" [0.6cm] and 3/32" [0.2cm] diameter burs, mandrels and chuck accessories with the slight movement of a handle while the motor is turning.

The Chuk/Changer provides the technician with exacting concentricity and quick, easy changing of various chucks which greatly speeds grinding, polishing and finishing procedures. The No.16 Chuk/Changer is adaptable to most types of dental motors, including older models e.g. **RED WING**, R & M, Ritter, Baldor, Buffalo, & Emesco. The 16 drastically reduces time you'll spend changing chucks and pays for itself quickly



16B CHUK/CHANGER

The No.16 is shipped complete with chuck accessories as follows. . .

P-16-15 Arbor Band Chuck	Qty-1
P-16-16 Spiral Tapered Chuck	Qty-1
P-16-19 Stone Chucks	Qty-2
P-16-20 Mandrels for small abrasive wheels.....	Qty-2
P-16-8 Collet Wrench	Qty-1
P-16-10 Clutch	Qty-1
P16-ADW Adapter	Qty-1

NOTE: Specify make, model and catalogue number of lathe when ordering. 16L Chuk/Changer is for use on left side of motor.



**MADE IN
THE USA**

HANDLER

OVER EIGHTY YEARS OF QUALITY MANUFACTURING AND SERVICE...

MANUFACTURING CO., INC. 612 NORTH AVENUE EAST, WESTFIELD, N.J. 07090

Ph: (908) 233-7796 Fax: (908) 233-7340 • 800-27-HANDL • E-Mail: info@handlermfg.com • web site: www.handlermfg.com

Sold Through Distributors Only

