

HANDLER

RED WING INTERNATIONAL

800.274.2635
www.handlermfg.com



Customer Testimonial

"My Handler 1/3 HP model trimmer is so strong, I used it to edge my outside walkway brick pavers and I still use it today in the lab. It is great."



For decades, Handler has been the trusted source in the dental industry for Red Wing Model Trimmers and now we have engineered them even better.

- **NEW & IMPROVED** – All Handler wet model trimmers now come standard with a pair of rubber suction cup feet and a reversible, REINFORCED Carbo-Dent wheel.
- The rubber suction cup feet lift the front for better water flow and add stability to the unit for improved safety.
- The 1/4" thick reinforced wheel for improved safety is stronger than most other wheels and gives you better peace of mind.
- All Handler trimmers are powered by high torque, long lasting maintenance free motors and trimmers come supplied ready for use.

Other standard features include:

- Easy access to the wheel.
Front doors swings open with the turn of a single thumb screw
- Full radius perforated, no splash water spray
- Adjustable, calibrated work table made of aluminum
- Aesthetically pleasing oven cured white powder coated finish
- Covered motor for improved safety and quieter operation
- Top mounted work tray to store cases for convenience when working

Models:

31X.....10" 1/4 HP

31.....10" 1/3 HP

32X.....12" 1/3 HP

32.....12" 1/2 HP

32-3/4...12" 3/4 HP

WE ARE ALL ABOUT QUALITY - WE ARE ALL ABOUT CHOICES

Available options and accessories ON ALL HANDLER TRIMMERS



#32SHLD

HANDLER designed the 32SHLD to provide protection against flying debris. This is a must for OSHA compliance. It has clear acrylic 6"x8"x1/8" shield and is attached to a magnetic base and a 13" adjustable arm for maximum flexibility.

#32-14INS

Factory installed Ortho Angularator precision machined Orthodontic Table

#32-15

Water solenoid valve for automatic water flow when trimmer is activated. It can be factory installed or installed in the field by a licensed plumber



#31KD / #32KD

12" diamond wheels. Various grit
10" & 12" Carbo-Dent reversible & reinforced 1/4" wheels any trimmer



31-DTS
DRY MODEL TRIMMER
Call for details

WWW.HANDLERMFG.COM

612 North Avenue East, Westfield N.J. 07090 Tel (908) 233-7796 Fax (908) 233-7340

Made in the U.S.A.
Sold Through Distributors Only