

# HANDLER

The Handler 31-DTS dry model trimmer system is designed to allow the user to trim plaster and stone models without water. The 3 stage system effectively and precisely trims models, while removing the dust from the work area of the model trimmer through a 5 gallon cyclone containment system.

The suction is provided by the powerful and popular Handler 62-II Super Sucker II. The Handler 31 model trimmer powered by a 1/3 HP sealed ball bearing, capacitor start motor with a 10" abrasive disc.



## 31-DTS DRY TRIMMER SYSTEM

### FEATURES:

- Handler 31 model trimmer with 10" abrasive disc, 5 gallon cyclone to trap dust particles for easy disposal
- A 62-II super sucker draws a powerful 110 cfm through the cyclone and entraps particles .5 microns and larger
- Hoses connect the 3 components, dust is removed from the work area & falls into a 5 gallon cyclone canister

### BENEFITS:

- Easy to set up and use
- Does not require water
- Any Handler model trimmer can easily be modified, components are interchangeable
- No more water hook ups, no more wet work areas, no more plaster traps, no more mess!
- Just empty the cyclone canister periodically in the normal trash & clean the cyclone cartridge.  
Change the 62-II filter bags when dirty
- Best of all any Handler model trimmer can be used wet or dry



**MADE IN  
THE USA**

CONTACT HANDLER MFG FOR MORE INFORMATION & HOW TO ORDER  
YOUR DRY SYSTEM TODAY!!!

**HANDLER**

**MANUFACTURING CO., INC.**

**OVER EIGHTY YEARS OF QUALITY MANUFACTURING AND SERVICE...**

612 NORTH AVENUE EAST, WESTFIELD, N.J. 07090

Ph: (908) 233-7796 Fax: (908) 233-7340 • 800-27-HANDL • E-Mail: [info@handlermfg.com](mailto:info@handlermfg.com) • web site: [www.handlermfg.com](http://www.handlermfg.com)

Sold Through Distributors Only

