

HANDLER

SUPER SUCKER II

#62II

110 CUBIC FEET/MIN.

QUICK CHANGE FILTER

1 YEAR WARRANTY

CONTINUOUS DUTY

SINGLE POSITION UNIT

ACTIVATED CHARCOAL FILTER

REDUCES MONOMER ODORS

THERMALLY PROTECTED MOTOR

QUIET/COMPACT DEVELOPS PEAK 1.4 HP



#62II

Developed to provide excellent suction for a single bench top work position. The 62-II Super Sucker develops 110cfm suction capacity at a static pressure of 63" of H₂O. Designed for continuous duty, its filter bag is easily changed through the front removable door. An activated charcoal filter helps remove monomer odors. Shipped with 2 1/2" flexhose and fishmouth.

Dimensions: 10" W x 15 1/4" D x 14" H

Sold Through Distributors Only

HANDLER

MANUFACTURING CO., INC.

OVER EIGHTY YEARS OF QUALITY MANUFACTURING AND SERVICE...

612 NORTH AVENUE EAST, WESTFIELD, N.J. 07090

Ph: (908) 233-7796 Fax: (908) 233-7340 • 800-27-HANDL • E-Mail: info@handlermfg.com • web site: www.handlermfg.com

