

86D SPLASH HOOD

**Model:**

86D Splash Hood

Description:

- New and improved sleek design
- Sturdy steel welded construction
- Attractive eagle white textured powder coated finish
- Steel insert pan - slides out for easy cleaning
- Bright, long lasting LED light with side mounted on/off switch
- Large, adjustable protective safety shield prevents splashing
- Rubber padded bottom for safety to prevent slipping and movement during procedures
- Comes fully assembled

Shield Dimensions: 6½" W x 7¾" H

Unit Dimensions: 6¼" W x 12" D x 12"

Box Dimensions: 14" W x 14" L x 9" H

Shipping weight: 9 pounds



Steel insert pan that slides in and out for easy cleaning



Long lasting, bright LED light



Rubber padded bottom to prevent slipping



Adjustable protective safety shield 6 ½"W x 7 ¾"H to prevent splashing on the user



Side mounted on/off switch for better illumination

Let HANDLER design and outfit your lab for efficiency and profitability. We can customize to fit any space.



THE RECOGNIZED
LAB PROFESSIONALS



908.233.7796 • HANDLERMFG.COM • MADE IN THE U.S.A.