

# HANDLER



## SYNCRONTORQUE IV 500-IV

# NEW



- ⚡ 713.8g/cm of torque provides ample power at all speeds
- ⚡ Up to 50,000 RPM for all applications including cutting, grinding, sanding, trimming, etching, carving, drilling and polishing
- ⚡ Brushless handpiece motor never requires brush changes, and offers minimal maintenance if cared for properly
- ⚡ Comes with a 1 year warranty against manufacturing defect
- ⚡ Reverse direction allows use with left or right hand accessories
- ⚡ Auto cruise lets you maintain any speed
- ⚡ Light weight and compact design with handle for portability
- ⚡ Foot pedal can be used on the floor or mounted as a knee treadle
- ⚡ Stylish & contemporary design for the modern lab
- ⚡ Affordable cost and quality construction for high value & optimum performance
- ⚡ Available in 230V, model 500-IVE

Shipping weight: 10 lbs.

Control box dimensions: 4 3/4" W x 4 3/4" H x 7 1/2" L

# HANDLER

[www.handlermfg.com](http://www.handlermfg.com)

MADE IN THE U.S.A

612 North Ave East  
PO Box 520  
Westfield, NJ 07090-0520  
(908) 233-7796  
(800) 274-2635