




HANDLER

RED WING INTERNATIONAL

575

Multifunctional Bench-Top Polishing Unit



575 Multifunctional Bench Top Polishing System

- Heavy gauge, stainless steel, powder coated cabinet
- Affordably priced significantly lower than competitive units.
- Easy operation for light, dry grinding, or polishing with pumice using the pumice pans
- Powerful suction system protects user with safer, healthier environment. (Helps to meet OSHA standards).
- Includes a 26A Red Wing Polishing lathe; the industry leader
- Easy access to reusable, replaceable pleated rear cartridge filters
- Replaceable, adjustable clear acrylic shield protects user from debris, and can be set at different heights to accommodate user sitting down or standing up
- Arm rests for comfortable, sustainable operation
- Designed & Built by HANDLER in the U.S.A.

115V/60Hz. Also available in 230V/50Hz and 230V/60Hz at additional cost.

Shipping Wgt: 115 lbs (Shipped in two cartons)

Unit Dimensions: 30 1/2" W x 16" H x 21" D

Also available without lathe - order Model 575-0

Shipping Wgt: 75 lbs.



575-0 Polishing Unit without Lathe

OVER 90 YEARS OF AMERICAN-MADE
DURABILITY IN DENTURE PROCESSING EQUIPMENT

WWW.HANDLERMFG.COM

612 North Avenue East, Westfield N.J. 07090 Tel (908) 233-7796 Fax (908) 233-7340

800.274.2635



Made in the U.S.A. 
Sold Through Distributors Only