

# HANDLER

RED WING INTERNATIONAL



*For decades, Handler has been the trusted source in the dental equipment industry.*

## **26105SSWC Stainless Steel Denture Processing System**

complete with vented cabinet allows the technician to Boil out, Wash out and Cure dentures. Designed and manufactured by Handler, the dual tank production system with separate temperature controls gives the technician the ability to cure in one tank while using the other tank to wash out or boil out. This processing method will increase efficiency in your lab and provide you with years of service.

**Don't trust your denture processing to anything less.**

## **26105SSWC** Boilout/Washout/Curing



- Dual tank system with separate operating controls
- Quality stainless steel tanks and top are seamlessly integrated
- Dual tank system is insulated for safety and to conserve energy
- Curing tank holds up to 16 flasks simultaneously
- State of the art electric heating elements with built in safety shut-off
- Heavy-duty pump and 4 faucet wax wash-out system
- Water-proof on/off switch to activate pump
- Grounded 6 feet - 110V cordset for pump
- Vented, heavy gauge steel cabinet with oven baked powder coated finish
- Water spill reducing backsplash and raised marine edge for work area safety
- Tank covers are included to bring water up to temperature more quickly
- Stainless steel wash out tray, flasks trays and one set of wax filters included
- Stand alone device for easy install with manual fill and manual drain
- Fill faucet (part #P26105-30F) available and can be easily mounted on the unit
- 26105SSWCCEL available in 208, 220 and 240 Volt, please specify when ordering
- 26105SSWCG is available in natural or propane gas, please specify when ordering
- Single tank, triple tank and custom tank systems available upon request
- Built to order and manufactured in the USA by Handler

### Available Colors:



**OVER 90 YEARS OF AMERICAN-MADE  
DURABILITY IN DENTURE PROCESSING EQUIPMENT**

**WWW.HANDLERMFG.COM** 800.274.2635  
612 North Avenue East, Westfield N.J. 07090 Tel (908) 233-7796 Fax (908) 233-7340



Made in the U.S.A.   
**Sold Through Distributors Only**