


NEW! PTL-1 PRO TECH LITE WORK STATION




Down Draft
Dust Collection
Portal with Shield



Deep Draws
for storage



Pop-up duplex
power supply with
2 USB ports



Stainless Steel
Surface Protector



PTL-1 Work Station for any Lab Application

"The 1 is a doer, a powerful force that produces results and does not allow anything or anyone to limit its potential. The 1 is aggressive, a necessary energy for creating and producing. The 1 is always in the forefront: a spear point directing and leading others."

- Modern practical design for any lab application
- 42" Wide x 24 1/2" Deep effective work surface
- Resilient all steel construction and solid work surface to meet demanding lab requirements
- Powerful, hidden and quiet 550A dust collector with disposable HEPA filter, easy to access and easy to clean
 - Available without 550A dust collector
 - Hidden suction inlet with screen and shield located on the work surface
- No need for armrest or attachments
- Two storage drawers (16" L x 13 1/2" W x 5 1/2" H) for tools and supplies
- Two tier shelf for case pans with retractable and swivel LED light
- Work station available without shelves and LED light
- Pop up duplex power supply with 2 USB ports
- Ample leg room for comfort
- Stainless Steel work surface protector
- Comes fully assembled, no installation required
- Plug into standard wall outlet
- Standard colors apply



550A Dust Collector Storage



NEW!

500-IV Syncro-Torque

Brushless, 50,000 RPM, Smooth, Hi-Torque, Forward and Reverse Hand Piece, Variable Speed Hand and Foot Control. Compact, Easy To Use, Affordable.