



HANDLER

REDWING INTERNATIONAL

PRO TECH *Lite*



PTL-231

Petite Plaster Bench

Shipping weight: 185lbs

Dimensions: 35" W x 36" H x 24" D

- Sturdy steel construction
- Leg levelers
- Durable oven-cured powder coated finish
- Light weight – 35"W x 24"D work surface
- Resistop® sink and drain
- Gooseneck lab faucet
- Optional in drawer plaster bin liner
- Optional kick plate hands free water activation
- Resistop® work surface
- Compact space design
- No assembly – easy install
- Affordable price

PTL-211 DESK

Shipping weight: 185lbs

Dimensions: 35" W x 36" H x 24" D



Available Colors:

MADE IN THE
U.S.A.
SINCE 1920

Warm Gray, Colonial White, Designer Beige,
Eagle White, Blue Heaven, Charcoal, Merlot

To order the PTL-231, please contact
your dealer representative.