

TECHNICAL DATA SHEET



PRODUCT

Ambi-ance® Ambidextrous Nitrile Gloves, Powder free

PRODUCT OVERVIEW

Ambidextrous Nitrile Gloves Powder-Free. Gloves are intended for use when conducting medical examinations, diagnostic and therapeutic procedures to protect patient and user from cross contamination/infection.

MATERIALS

Nitrile (Carboxylate Butadiene Acrylonitrile)

DETAILS

- Maximum comfort and performance with textured fingertips for improved dexterity, enhanced grip and tactile sensitivity.
- Soft stretch Nitrile for a comfortable custom fit.
- Puncture and tear resistant barrier.
- Beaded cuff for extra protection.

PRODUCT/ PACKAGING

Product is available in the following configurations

Product Code	Description
19010	Ambi-ance® Ambidextrous Nitrile Gloves, Extra Small, 200/box
19011	Ambi-ance® Ambidextrous Nitrile Gloves, Small, 200/box
19012	Ambi-ance® Ambidextrous Nitrile Gloves, Medium, 200/box
19013	Ambi-ance® Ambidextrous Nitrile Gloves, Large, 200/box
19014	Ambi-ance® Ambidextrous Nitrile Gloves, Extra Large, 200/box

PRODUCT SPECIFICATIONS

SIZE	X-SMALL	SMALL	MEDIUM	LARGE	X-LARGE
Length (mm)	235	235	235	235	235
Width (mm)	75±5	85±5	95±5	106±5	115±5
Palm Thickness (mm) *Single Wall	0.06-0.09	0.06-0.09	0.06-0.09	0.06-0.09	0.06-0.09
Finger (mm)	0.08-0.14	0.08-0.14	0.08-0.14	0.08-0.14	0.08-0.14
Beading Thickness	≥0.5	≥0.5	≥0.5	≥0.5	≥0.5
Cuff Thickness (mm)	0.05-0.08	0.05-0.08	0.05-0.08	0.05-0.08	0.05-0.08

STORAGE REQUIREMENTS

Store in a dry cool place; avoid Ozone.

SHELF LIFE

Three (3) years from the date of manufacture, when above storage conditions are followed.

Crosstex® and Ambi-ance® are registered trademarks of Crosstex International, Inc.