

Crosstex is now a proud member of HuFriedyGroup



Technical Data Sheet

Dental Dams



| Product | Latex – Powder Free/Low Protein Dental Dams | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|------------------------|---|--|--|--------------|-------------|----------------------|--|-------|--|-------|--|-------|---|-------|--|-------|--|-------|---|------------------|--|-------|--|-------|--|-------|---|-------|--|-------|--|-------|---|
| Product Overview | Used to isolate and protect teeth and mouth during a dental procedure | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Materials | Manufactured from medium modulus natural latex. | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Details | <ul style="list-style-type: none">• High tear-strength• 100% Non-chlorinated• Powder-free with less than 5mg/sheet.• Low protein: extractable protein content is below FDA allowable limit of 50µg/sheet (ASTM D5712) | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Storage Requirements | Keep dry and away from sunlight | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Shelf Life | 5 years | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Packaging | <table><tr><th>Product Code</th><th>Description</th></tr><tr><td colspan="2">Mint flavored, Green</td></tr><tr><td>19202</td><td>5 in x 5 in, Thin gauge, 52 Sheets/Box</td></tr><tr><td>19200</td><td>5 in x 5 in, Medium gauge, 52 Sheets/Box</td></tr><tr><td>19201</td><td>5 in x 5 in, Heavy gauge, 52 Sheets/Box</td></tr><tr><td>19402</td><td>6 in x 6 in, Thin gauge, 36 Sheets/Box</td></tr><tr><td>19400</td><td>6 in x 6 in, Medium gauge, 52 Sheets/Box</td></tr><tr><td>19401</td><td>6 in x 6 in, Heavy gauge, 52 Sheets/Box</td></tr><tr><td colspan="2">Unflavored, Blue</td></tr><tr><td>19102</td><td>5 in x 5 in, Thin gauge, 52 Sheets/Box</td></tr><tr><td>19100</td><td>5 in x 5 in, Medium gauge, 52 Sheets/Box</td></tr><tr><td>19101</td><td>5 in x 5 in, Heavy gauge, 52 Sheets/Box</td></tr><tr><td>19302</td><td>6 in x 6 in, Thin gauge, 36 Sheets/Box</td></tr><tr><td>19300</td><td>6 in x 6 in, Medium gauge, 52 Sheets/Box</td></tr><tr><td>19301</td><td>6 in x 6 in, Heavy gauge, 52 Sheets/Box</td></tr></table> | | | Product Code | Description | Mint flavored, Green | | 19202 | 5 in x 5 in, Thin gauge, 52 Sheets/Box | 19200 | 5 in x 5 in, Medium gauge, 52 Sheets/Box | 19201 | 5 in x 5 in, Heavy gauge, 52 Sheets/Box | 19402 | 6 in x 6 in, Thin gauge, 36 Sheets/Box | 19400 | 6 in x 6 in, Medium gauge, 52 Sheets/Box | 19401 | 6 in x 6 in, Heavy gauge, 52 Sheets/Box | Unflavored, Blue | | 19102 | 5 in x 5 in, Thin gauge, 52 Sheets/Box | 19100 | 5 in x 5 in, Medium gauge, 52 Sheets/Box | 19101 | 5 in x 5 in, Heavy gauge, 52 Sheets/Box | 19302 | 6 in x 6 in, Thin gauge, 36 Sheets/Box | 19300 | 6 in x 6 in, Medium gauge, 52 Sheets/Box | 19301 | 6 in x 6 in, Heavy gauge, 52 Sheets/Box |
| Product Code | Description | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Mint flavored, Green | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 19202 | 5 in x 5 in, Thin gauge, 52 Sheets/Box | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 19200 | 5 in x 5 in, Medium gauge, 52 Sheets/Box | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 19201 | 5 in x 5 in, Heavy gauge, 52 Sheets/Box | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 19402 | 6 in x 6 in, Thin gauge, 36 Sheets/Box | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 19400 | 6 in x 6 in, Medium gauge, 52 Sheets/Box | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 19401 | 6 in x 6 in, Heavy gauge, 52 Sheets/Box | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Unflavored, Blue | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 19102 | 5 in x 5 in, Thin gauge, 52 Sheets/Box | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 19100 | 5 in x 5 in, Medium gauge, 52 Sheets/Box | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 19101 | 5 in x 5 in, Heavy gauge, 52 Sheets/Box | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 19302 | 6 in x 6 in, Thin gauge, 36 Sheets/Box | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 19300 | 6 in x 6 in, Medium gauge, 52 Sheets/Box | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 19301 | 6 in x 6 in, Heavy gauge, 52 Sheets/Box | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Regulatory Information | Product Code: EIE Classification: Class 1, 510(k) Exempt Regulation: 872.6300, Rubber Dam and Accessories Medical Device Listing: D142786 Single-use | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

All product names are trademarks of Crosstex International, Inc., a Cantel Medical Company, its affiliates or related companies, unless otherwise noted.
 © 2020 Crosstex International, Inc. DLIT01576 Rev A 0820

Crosstex International, Inc.

Phone 888-276-7783 | Fax 631-582-1726

10 Ranick Road | Hauppauge, NY 11788-4209