

TECHNICAL DATA SHEET



PRODUCT Sparkle™ Disposable Air Water Syringe Tips

PRODUCT OVERVIEW The use of forced air and water in dentistry is an everyday occurrence. During exams, prophylaxis/scaling appointment, restorative care, etc., it is necessary to flush or rinse debris away from a tooth or dry the tooth with forced air. This is accomplished using the delivery system and the air water syringe. The syringe needs a tip to deliver the air and water as required.

MATERIALS The Air Water Syringe Tips are made of PVC and Tritan MX811.

DETAILS

- Quad Core Air Channel ensures air/water separation.
- Beveled edge provides increased patient comfort.
- Single-use eliminates cross-contamination risks and the cleaning process required for reusable tip.
- Tip connects securely without the need for an adapter.
- Rigid construction provides safe and reliable cheek retraction.
- Translucent material allows visualization of water droplets, flow and/or blockages.

PRODUCT/ Product is available in the following configurations

PACKAGING

Product Code	Description
BCSAWSCL	Clear, 250/Bag
BCSAWSWH	White, 250/Bag
BCSAWS	Assorted (Blue, Purple, Yellow, Amber), 250/Bag
BCSAWSB	Assorted (Blue, Purple, Yellow, Amber), 1500/Bag

STORAGE REQUIREMENTS The packaged products should be stored in an area with a relative humidity no more than 80%, no corrosive gas and good ventilation.

SHELF LIFE Five years

0818 DTDS00326

Crosstex™ and Sparkle™ are trademarks of Crosstex International, Inc., a Cintel Medical Company