

# INSTRUCTIONS FOR USE



## Operatory Protection

### Syringe Sleeve

#### Instructions for use:

##### With opening:

1. Plastic sleeve has opening at end.
2. Insert syringe into sleeve.

##### Without opening:

1. Place plastic sleeve over syringe and push through.



Ref.#'s BCSS (without opening),  
BCSSO, BCSSSBO (with opening)



Do not reuse.



Not made with natural rubber latex.



**CROSSTEX**<sup>®</sup>  
A CANTEL MEDICAL COMPANY



MDSS GmbH  
Schiffgraben 41  
30175 Hannover, Germany



Made in USA

Manufactured by: Crosstex International, Inc.  
10 Ranick Road | Hauppauge, NY 11788 | 631.582.6777 | 888.276.7783 | [Crosstex@crosstex.com](mailto:Crosstex@crosstex.com)

BCSS-IFU  
Rev. 0 – 5/17