

INSTRUCTIONS FOR USE



Operatory Protection

Syringe Sleeve with Adhesive Strip

Instructions for use:

1. Plastic sleeve has opening at end.
2. Insert syringe into sleeve.
3. Seal closed with adhesive to create a closed barrier.



Ref.#'s BCSSA



Do not reuse.



Not made with natural rubber latex.



CROSSTEX[®]
A CANTEL MEDICAL COMPANY



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



Made in China

Manufactured for: Crosstex International, Inc.
10 Ranick Road | Hauppauge, NY 11788 | 631.582.6777 | 888.276.7783 | Crosstex@crosstex.com

BCSSA-IFU
Rev. 0 – 5/17