

INSTRUCTIONS FOR USE



Operatory Protection

Tray Sleeves

Instructions for use:

1. Insert tray into plastic sleeve.
2. Fold over plastic end for a secure fit.



Ref.#'s BCTS1014, BCTS710, BCTS1114, BCTS1116



Do not reuse.



Not made with natural rubber latex.



CROSSTEX
A CANTEL MEDICAL COMPANY



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



Made in USA

Manufactured by: Crosstex International, Inc.
10 Ranick Road | Hauppauge, NY 11788 | 631.582.6777 | 888.276.7783 | Crosstex@crosstex.com

BCT-IFU
Rev. 0 – 5/17