

INSTRUCTIONS FOR USE

Imaging Comfort




Sensor Slippers® Digital X-Ray Comfort Barriers

How to use your Crosstex® Sensor Slippers®:



Step 1:
Insert digital sensor into
Sensor Slipper.



Step 2:
Insert cushioned digital sensor
into gripping-style holder of
your choice.



Step 3:
Take images and discard Sensor
Slipper after use.

Ref.#'s BWSSL

Sensor Slippers® is a registered trademark of Crosstex International, Inc.



Not made with natural rubber latex.



CROSSTEX®
A CANTEL MEDICAL COMPANY

Manufactured for: Crosstex International, Inc.

10 Ranick Road | Hauppauge, NY 11788 | 631.582.6777 | 888.276.7783 | Crosstex@crosstex.com

Made in China

BWSS-IFU
Rev. 0 – 11/16