

INSTRUCTIONS FOR USE

Treatment Accessories



Medicine/Mixing Cups

Instructions for use:

- Convenient for dispensing or mixing.
- Features ounce, cc, ml, dram & teaspoon calibrations.
- Unbreakable plastic; can be squeezed to form pour spout.
- Multi-purpose product.



Ref.#'s CX1



Do not reuse.



Not made with natural rubber latex.



CROSSTEX[®]
A CANTEL MEDICAL COMPANY



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



Made in USA

Manufactured for: Crosstex International, Inc.
10 Ranick Road | Hauppauge, NY 11788 | 631.582.6777 | 888.276.7783 | Crosstex@crosstex.com

CX1-IFU
Rev. 0 – 11/16