

Crosstex™ Disposable Face Shields

Reduce your risk of exposure by protecting your critical zone.

- **Anti-fog** treatment on both the inside and outside of the shield for optimum visibility.
- High-quality, **made-in-USA** film offers maximum optical clarity.
- 1.5 inch vented foam headband holds shield away from face to **improve airflow** while exerting no pressure on temples.
- Special ultrasonically welded elastic headband for **added strength**.
- Smooth, beveled edges for **increased safety**.
- **Single-use** minimizes cross-contamination.
- **Individually wrapped** to protect the integrity and cleanliness of each shield.



Pictured
¾ length 7 in
Clear Face Shield



Use Crosstex™ Disposable Face Shields during procedures likely to generate splashes, spatter or aerosol

The central areas of the face are the most vulnerable for transmission especially during certain procedures. Use layered PPE like glasses, a mask and a face shield to increase the protection in these areas.

The Crosstex™ Disposable Face Shield protects your mask and face from direct spatter and may prolong mask life.¹

CAUTION: Outside of masks and face shields are likely to become contaminated during use. Wash hands after touching any contaminated surfaces.

- ✓ CDC recommends that a face shield be included as part of the PPE for dental healthcare personnel during all procedure types to minimize risk of COVID-19 transmission.²
- ✓ OSHA requires that all employers supply eyewear or face shield protection for all employees.
- ✓ AAMI ST79 recommends the use of face shields during manual cleaning.

ORDERING INFORMATION

Ref. #	Description	Quantity
GCSS	¾ length 7 in Clear Face Shield	24/Box
GCSSB	¾ length 7 in Clear Face Shield, bulk	72/Box
GCSS9	Full length 9 in Clear Face Shield	18/Box
GCSS9B	Full length 9 in Clear Face Shield, bulk	72/Box

Not made with natural rubber latex.

1) Face shields do not replace the use of a procedure-appropriate face mask. Do not use a face shield without a face mask during aerosol procedures. 2) CDC Guidance for Dental Settings (Aug. 28, 2020)
All company and product names are trademarks of Crosstex International, Inc., its affiliates or related companies, unless otherwise noted.
©2020 Crosstex International, Inc. All rights reserved. DLIT00081 Rev B 1020