

TECHNICAL DATA SHEET



PRODUCT Crosstex® Articulating Film

PRODUCT OVERVIEW Marks high and low occlusion points and provides the most accurate way to mark every point of articulation. Resulting in clear positive marking without smudging on dry or wet fields.

- MATERIALS**
- Plastic Film with wax pigmented coating.
 - Not made with natural rubber latex.

- DETAILS**
- Superior manufacturing process combined with high quality control guarantees clear, positive marking without smudging on wet or dry fields.
 - Color coded packaging and product descriptions make product identification easy.
 - Articulating film available in blue, red, blue/red combo and black/red combo.
 - Available in single sided and double sided straight formats.

**PRODUCT/
PACKAGING**

Product is available in the following configurations:

Product Code	Description
	All packages contain 25 sheets per book and 12 books per box
TFMBKR	Black/Red Combo, Precisely True, 0.0005"/ 13 µm
TPP	Blue, Thin Precisely True, 0.0005"/ 13 µm
TFMR	Red Precisely True, 0.0005"/ 13 µm
TFMBR	Red/Blue Combo, Precisely True, 0.0005"/ 13 µm

STORAGE REQUIREMENTS Store in dry, room temperature location.

SHELF LIFE Single use only. Professional use only.

MKT421
Crosstex® is a registered trademark of Crosstex International, Inc