

TECHNICAL DATA SHEET



PRODUCT

Crosstex® Articulating Foil

PRODUCT OVERVIEW

Marks high and low occlusion points and provides the most accurate way to mark every point of articulation. Resulting in clear positive marking without smudging on dry or wet fields.

MATERIALS

- Mirrored Plastic Film with wax pigmented coating.
- Made with Polyethylene Terephthalate (PET).

DETAILS

- Superior manufacturing process combined with high quality control guarantees clear, positive marking without smudging on wet or dry fields.
- Color coded packaging and product descriptions make product identification easy.
- Articulating foil available in red, blue and red/blue combination.
- Available in single sided and double sided straight formats.

PRODUCT/ PACKAGING

Product is available in the following configurations:

Product Code	Description
	All packages contain 25 sheets per book and 12 books per box
TFOBR	Red/Blue Combo, 0.0005"/12 µm
TFOR2	Red Two-Sided, 0.0005"/12 µm
TFOB2	Blue Two-Sided, 0.0005"/12 µm
TFOB	Blue One-Sided, 0.0005"/12 µm

STORAGE REQUIREMENTS

Store in dry, room temperature location.

SHELF LIFE

Single use only. Professional use only.

MKT422
Crosstex® is a registered trademark of Crosstex International, Inc.