

Safety Data Sheet

VioNexus

Section 1

Product Description

Product Name:	VioNexus Antimicrobial Foam Soap
Recommended Uses:	Hand Hygiene
Synonyms:	Foaming Antimicrobial Hand Soap
Distributor:	Metrex Research
	28210 Wick Road Romulus, MI 48174
Chemical Information	
Emergency:	
Chemtrec	1.800.424.9300

Section 2

Hazard Identification

OSHA Regulatory Status--This product contains no hazardous components as defined in the 2012 OSHA Hazard Communication, 29 CFR §1910.1200. Warnings below are for exposure to large quantities of the product. Normal usage should not create hazardous conditions.

Appearance—Aqueous solution

Physical state—Liquid

Odor—Mild Odor

Can cause eye irritation.

Precautionary Statements—Prevention

Keep away from heat/sparks/open flames/hot surfaces-----No Smoking.

Use only non-sparking tools.

Keep container tightly closed.

Precautionary Statements—Response

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.

IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin (hair) with water/shower.

IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.

In case of fire: Use CO₂, dry chemical, or foam for extinction.

Precautionary Statements—Storage

Stored in a well-ventilated place. Keep Cool.

Safety Data Sheet

Precautionary Statements—Disposal

Dispose of contents/container to an approved waste disposal plant.

Hazards not otherwise classified (HNOC)

Not Applicable

Other Information

Unknown Acute Toxicity

Section 3

Composition/Information on Ingredients

Chemical Name	CASNo.	Weight-%
4-Chloro-3, 5 xylene (Chloroxylene)	88-04-0	< 2

Section 4

First Aid Measures

First Aid Measures

Eye Contact

Hold eye(s) open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after first 5 minutes, then continue rinsing eye(s). If eye irritation persists: Get medical advice/attention.

Skin Contact

For accidental exposure to large quantities, wash off immediately with soap and plenty of water while removing all contaminated clothes and shoes. Wash contaminated clothing before reused. Get medical attention if irritation develops.

Inhalation

Remove to fresh air. If signs/symptoms continue, get medical attention.

Ingestion

Do NOT induce vomiting. Never give anything by mouth to an unconscious person. Call physician.

Most important symptoms and effects, both acute and delayed

Symptoms

See Section 11 for symptom information.

Indication of any immediate medical attention and special treatment needed

Note to physicians

Treat symptomatically.

Section 5

Fire-Fighting Measures

Suitable Extinguishing Media

Dry Chemical, Water spray (fog), Carbon dioxide (CO₂), Foam.

Unsuitable Extinguishing Media

No Information available.

Specific hazards arising from the chemical

No Information available.

Hazardous combustion products

Carbon Monoxide. Carbon Dioxide (CO₂).

Explosion Data

Sensitivity to Mechanical Impact

None

Sensitivity to Static Discharge

May be ignited by friction, heat, sparks or flames.

Protective Equipment and Precautions for Firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective

Safety Data Sheet

gear. Cool containers with flooding quantities of water until well after fire is out.

Section 6

Accidental Release Measures

Personal Precautions, Protective Equipment and Emergency Procedures

Personal Precautions Use personal protection recommended in Section 8. Ensure adequate ventilation, especially in confined areas.

For emergency responders Isolate area. Keep unnecessary personnel away.

Environment Precautions

Environmental Precautions Prevent entry into waterways, sewers, basements or confined areas. See Section 12 for additional ecological information.

Methods and material for containment and cleaning up

Methods for containment Prevent further leakage or spillage if safe to do so. Contain and collect spillage with non-combustible absorbent material, (e.g. sand, earth, diatomaceous earth, vermiculite) and place in container for disposal according to local/national regulations (See Section 13).

Methods for cleaning up Use clean non-sparking tools to collect absorbed material. May be ignited by friction, heat, sparks or flames. Collect spillage. Soak up with inert absorbent material. Sweep up and shovel into suitable containers for disposal. Following product recovery, flush area with water.

Section 7

Handling and Storage

Precautions for Safe Handling

Advice on Safe Handling Keep away from heat, sparks, flame and other sources of ignition (i.e., pilot lights, electric motors and static electricity). Take precautionary measures against static discharges. Do not eat, drink or smoke when using this product. Use personal protection recommended in Section 8. Avoid contact with skin, eyes or clothing. Use only in well-ventilated areas. Avoid breathing vapors or mists. Wash thoroughly after handling. Handle in accordance with good industrial hygiene and safety practice.

Conditions for safe storage, including any incompatibilities

Storage Conditions Keep Containers tightly closed in a dry, cool and well-ventilated place. Keep from freezing. Keep away from heat, sparks, flame and other sources of ignition (i.e., pilot lights, electric motors and static electricity). Store between 40 to 100°F.

Incompatible materials Heat, sparks, open flame, other ignition sources. Reacts violently with strong oxidants such as nitric acid and silver nitrate causing fire and explosion hazard. Reacts slowly with calcium hypochlorite and ammonia causing fire and explosion hazard.

Safety Data Sheet

Section 8

Protection Information

Exposure Guidelines

Appropriate Engineering Controls

Engineering Controls Showers, eyewash stations, ventilation system.

Individual Protection Measures, such as personal protective equipment

Eye/Face protection	Wear safety glasses with side shields (or goggles).
Skin and body protection	Wear protective Neoprene™ gloves. Rubber gloves. Normal work clothing (long sleeved shirt and long pants) is recommended).
Respiratory protection	If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be provided in accordance with current local regulations.
General Hygiene Considerations	Wash face, hands and any exposed skin thoroughly after handling. Wash contaminated clothing and shoes before reuse. Do not Eat, Drink or Smoke when using this product.

Section 9

Physical and Chemical Properties

Information on basic physical and chemical properties

Formula:	See Section 3	Physical State:	Liquid
Odor:	Mild	Appearance:	Aqueous solution
Odor Threshold:	No Information Available	Color:	Green

<u>Property</u>	<u>Values</u>
pH	9.0-11.0
Melting Point/Freezing Point	No Information available
Boiling Point/ Boiling Range	212°F
Flash Point	No Information available
Evaporation rate	No information available
Flammability (solid, gas)	No information available

Safety Data Sheet

Flammability Limit in Air

Upper flammability limit:

No Information available

Lower flammability limit:

No Information available

Vapor pressure:

No Information available

Vapor density:

No information available

Specific Gravity

1.02

Water solubility

Soluble

Partition coefficient

No information available

Autoignition temperature

No information available

Decomposition temperature

No information available

Kinematic viscosity

No information available

Dynamic viscosity

No information available

Explosive properties

No information available

Oxidizing properties

No information available

Section 10

Stability and Reactivity Data

Reactivity

No data available.

Chemical Stability

Stable under recommended storage conditions.

Possibility of Hazardous Reactions

Vapors may form explosive mixture with air.

Conditions to avoid

Heat, flames and sparks, extreme temperature conditions

Incompatible materials

Strong oxidizing agents, strong acids, heat, sparks, open flame, other ignition sources.

Hazardous Decomposition Products

Carbon monoxide. Carbon dioxide (CO₂). Nitrogen.

Section 11

Toxicity Data

Information on likely routes of exposure

Product Information

Inhalation

No irritant under normal use.

Eye Contact

Can be irritating to eyes.

Skin Contact

No irritant under normal use.

Ingestion

May be harmful if swallowed.

Chemical Name	Oral LD50	Dermal LD50	Inhalation LC50
4-Chloro-3, 5 xyleneol (Chloroxylenol) 88-04-0	500 mg/kg	> 2000 mg/kg (Rat)	> 6.29 mg/L (Rat)

Safety Data Sheet

Information on toxicological effects

Symptoms No Information Available

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Sensitization No Information Available

Germ cell mutagenicity No Information Available

Carcinogenicity Not classified

Reproductive toxicity No Information Available

STOT single exposure No Information Available

STOT repeated exposure No Information Available

Aspiration hazard No Information Available

Section 12

Ecological Data

Ecotoxicity

Chemical Name	Algae/Aquatic Plants	Fish	Crustacea
4-Chloro-3, 5 xlenol (Chloroxylenol) 88-04-0	-	0.76 mg/L 96 h Oncorhynchus mykiss LC50	7.7 mg/L 48 h Daphnia magna mg/L EC50

Safety Data Sheet

Persistence and degradability

No Information Available.

Bioaccumulation

Mobility

No information available.

Other adverse effects

No information available

Section 13

Disposal Information

Waste treatment methods

Disposal of wastes

Disposal should be in accordance with applicable regional, national and local laws and regulations.

Contaminated packaging

Dispose of in accordance with federal, state and local regulations.

Section 14

Transport Information

DOT

Non-hazardous

Safety Data Sheet

Section 15

Regulatory Information

International Inventories

TSCA	Complies
DSL/NDSL	Complies
EINECS/ELINCS	Complies

Legend:

TSCA—United States Toxic Substances Control Act Section 8(b) Inventory

DSL/NDSL—Canadian Domestic Substances List/Non-Domestic Substances List

EINECS/ELINCS—European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances

US Federal Regulations

SARA 313

Section 313 of Title III of the Superfund Amendments and Reauthorization of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

SARA 311/312 Hazard Categories

Acute health hazard	No
Chronic Health Hazard	No
Fire Hazard	No
Sudden Release of pressure hazard	No
Reactive Hazard	No

CWA (Clean Water Act)

This product does not contain any substance regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42)

CERCLA

This material, as supplied, does not contain any substance regulated as a hazardous substance under the Comprehensive Environment Response Compensation and Liability Act (CERCLA) (40 CFR 302) or the Superfund Amendments and Reauthorization Act (SARA) (40 CFR 355). There may be specific reporting requirements at the local regional or state level pertaining to releases of this material.

US State Regulations

California Proposition 65

This product does not contain any Proposition 65 chemicals.

U.S. EPA Label Information

EPA Pesticide Registration Number	Not Applicable
--	----------------

Safety Data Sheet

Section 16

Additional Information

<u>HMIS</u>	Health Hazards	Flammability	Physical Hazards	Personal protection
	1	0	0	B (safety glasses,gloves)

Prepared by: Technical Department
Revision Date May 15,2015
Version 0
Revision Note No information available

Disclaimer

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

End of Safety Data Sheet

LS 104-368-00.000