



TheSharpChoice®

NMSDC Certified Diversity Supplier

DATASHEET

RELI® Sutures



MYCO Medical's Reli® brand sutures bring meaning to the definition of Reli·able: high quality at an affordable price.

- Knowledgeable team providing assistance for all your suturing needs
- REDIPASS™ needle-coating technology enables smoother needle penetration
- Complete range of sterile absorbable and non-absorbable suture products
- Available in a variety of suture sizes and lengths, and needle types and sizes
- Sutures packaged in an easy-retrieval pack that enables convenience and speed upon opening
- Every suture strand and needle attachment is tested to ensure a secure bond



Needles are ready to arm directly after opening the suture pack

Reli® Absorbable Sutures

All Reli® Sutures are available packaged as 12 packs/box, 30 boxes/case.

REDIGUT® PLAIN

NEEDLE				MATERIAL					
Name	Point	Curve	Length	Color	Length	6-0	5-0	4-0	3-0
MP-3	PPT REV CUT	3/8c	13mm	Natural	18"	—	GP686-BRC	—	—
MC-2	REV CUT		19mm		10"	—	—	GP531-BRC	—
MFS-2			24mm		18"	—	—	GP7771-BRC	GP8881-BRC
MFS-1			27"		—	GP820-BRC	GP821-BRC	GP822-BRC	
			18"		—	—	—	GP324-BRC	
MX-1	1/2c	23mm	27"	—	—	—	GP3247-BRC		
					18"	—	—	—	GP612-BRC

REDIGUT® CHROMIC

NEEDLE				MATERIAL					
Name	Point	Curve	Length	Color	Length	6-0	5-0	4-0	3-0
MP-3	PPT REV CUT	3/8c	13mm	Natural	18"	—	GC687-BRC	GC1654-BRC	—
MC-2	REV CUT		19mm			GC621-BRC	—	—	
MFS-2						—	—	GC6354-BRC	GC6364-BRC
MFS-1			24mm		—	GC634-BRC	GC635-BRC	GC636-BRC	
MX-1			1/2c		23mm	—	—	GC421-BRC	—

REDICAPRONE™ — 25 (PGCL)

for short/medium-term wound support; comparable to *Monocryl*

NEEDLE				MATERIAL						
Name	Point	Curve	Length	Construction & Color		Length	4-0	3-0	2-0	1
MFS-2	REV CUT	3/8c	19mm	Monofilament	Undyed	36"	MS922-SI	—	—	—
MFS-1			24mm				—	MS942-SI	MS943-SI	—
MSH	TAP PT	1/2c	26mm			27"	MS315-SI	MS316-SI	—	—
MCT-1			37mm				36"	—	—	MS732-SI
						—		—	—	MS347-SI

REDISORB™ (PGA)

for medium-term wound support; comparable to *Vicryl*

NEEDLE				MATERIAL						
Name	Point	Curve	Length	Construction & Color	Length	6-0	5-0	4-0	3-0	
MP-1	PPT REV CUT	3/8c	11mm	Braided	Violet	18"	V489-SI	—	—	
MP-3			13mm				27"	—	V463-SI	—
MPS-2					19mm			Undyed	VU492-SI	VU493-SI
MC-2			13mm				18"		—	—
MFS-2	REV CUT				19mm	Violet		—	V391-SI	V392-SI
MFS-1			24mm				Undyed	—	VU421-SI	VU422-SI
					Violet	18"		—	—	—
						Undyed	27"	—	—	VU441-SI

REDIDIOX™ (PDO)

for medium/long-term wound support; comparable to *PDS II*

NEEDLE				MATERIAL						
Name	Point	Curve	Length	Construction & Color		Length	4-0	3-0	2-0	0
MFS-2	REV CUT	3/8c	19mm	Monofilament	Violet	27"	ML397-SI	ML398-SI	—	—
MFS-1			24mm				—	ML452-SI	ML451-SI	—
MCP		1/2c	40mm				—	—	ML195-SI	—
MCP-2			26mm				—	—	—	ML970-SI
MCP-1			36mm				—	—	ML466-SI	—
MSH	TAP PT	1/2c	26mm				—	ML316-SI	—	—
MCT-2			26mm				—	—	ML333-SI	—



Reli® Non-Absorbable Sutures

All Reli® Sutures are available packaged as 12 packs/box, 30 boxes/case.

REDISILK™ (Silk) comparable to <i>Perma-Hand Silk</i>											
NEEDLE				MATERIAL							
Name	Point	Curve	Length	Construction & Color	Length	6-0	5-0	4-0	3-0	2-0	
MP-3	PPT REV CUT	3/8c	13mm	Braided	Black	18"	—	SK640-SI	SK641-SI	—	
MPS-3			16mm				SK681-SI	—	—	—	
MFS-2	19mm		—				SK682-SI	SK683-SI	SK7772-SI	—	
MFS-1	24mm		—				—	—	SK684-SI	—	
MFS	26mm		—				—	—	—	SK685-SI	
MJ-1	REV CUT	1/2c	16mm				—	—	—	SK7762-SI	—
MX-1			23mm				—	—	—	SK632-SI	—
REDILON™ (Nylon) comparable to <i>Ethilon</i>											
NEEDLE				MATERIAL							
Name	Point	Curve	Length	Construction & Color	Length	6-0	5-0	4-0	3-0	2-0	
MP-1	PPT REV CUT	3/8c	11mm	Monofilament	Black	18"	N697-SI	—	—	—	
MP-3			13mm				N1698-SI	N698-SI	N699-SI	—	
MPS-3			16mm				N1665-SI	N1668-SI	—	—	
MPS-2			19mm				—	N1666-SI	N1667-SI	—	
MPS-1			24mm				—	—	N1663-SI	—	
MC-2	REV CUT		13mm			10"	N667-SI	N6681-SI	N931-SI	—	
MFS-3			16mm				—	N668-SI	—	—	
MFS-2			19mm				10"	N660-SI	—	—	—
MFS-1			24mm					—	—	N6621-SI	—
MFS			26mm					—	N661-SI	N662-SI	N928-SI
			18"	—	—	—	N663-SI	—			
				—	—	—	—	—	N664-SI		
REDILENE™ (Polypropylene) comparable to <i>Prolene</i>											
NEEDLE				MATERIAL							
Name	Point	Curve	Length	Construction & Color	Length	6-0	5-0	4-0	3-0	2-0	
MP-1	PPT REV CUT	3/8c	11mm	Monofilament	Blue	18"	P8697-SP	—	—	—	
MP-3			13mm				P8695-SP	P8698-SP	P8699-SP	—	
MFS-3	REV CUT		16mm				P8660-SP	—	—	—	
MFS-2			19mm				—	P8661-SP	P8683-SP	—	
MFS-1			24mm				—	—	—	P8684-SP	
REDIBOND™ (Polyester) comparable to <i>Ethibond</i>/<i>Ti-Cron</i>											
NEEDLE				MATERIAL							
Name	Point	Curve	Length	Construction & Color	Length	6-0	5-0	4-0	3-0	2-0	
MC-2	REV CUT	3/8c	13mm	Braided	Green	18"	—	—	—	PT1001-SI	
MFS-2			19mm				—	—	—	PT1002-SI	
MV-5	TAP CUT	1/2c	18mm				—	PT7773-SI	—	—	



As a leading NMSDC Diversity Certified Supplier, MYCO Medical is committed to “Bringing Reliability, Innovation and Value” to the Healthcare industry with an emphasis on “physician preference” surgical devices and disposables.

To order, learn more or receive complimentary evaluation samples:

Please call MYCO Medical at (800) 454-6926 or (919) 460-2535 or email us at info@mycomedical.com

See the complete line of MYCO innovative disposable and reusable medical devices at www.mycomedical.com



MYCO
MEDICAL
TheSharpChoice®

2015 Production Drive
Apex, NC 27539 USA

info@mycomedical.com

www.mycomedical.com

Some medical devices are not licensed for sale in Canada.

© 2017-2019 MYCO Medical, Inc. All rights reserved.

LIT-00346 Rev 1 Eff 10.02.2019