



**Works Better.
Lasts Longer.
Costs Less.**



IMPRESSION



Bite Registration Materials Setting Times					
Viscosity	Color	Work Time	Set Time	Total Set Time	
Bite (unflavored)	Regular Set	2:00	1:00	3:00	
	Fast Set	Lt Blue :30	1:00	1:30	
	Super Fast Set	:30	:30	1:00	
Mint Choc Chip	Fast Set	Lt Blue :30	1:00	1:30	
	Super Fast Set	:30	:30	1:00	
Hard Bite (unflavored)	Fast Set	Blue :30	1:00	1:30	

“Works Better”

With excellent handling properties and accuracy this material provides the ideal viscosity and flow under pressure for precise bite registrations.

“Lasts Longer”

It stays where you put it – can be expressed directly onto the occlusal surface.

“Costs Less”

The superior formula saves on "re-takes" and is still about half the cost of the leading competitor.

Bite Registration Material

This Mousse-like and Thixotropic, fluffy vinyl polysiloxane material offers the best option to be used for bite registration, copying transfers, and primary impression with bite trays. The thixotropic material has an excellent durometer of over 80 and flows easily, yet will not run or sag.

Packaging:	Item No:	Compares to*:
2 cartridges (50 ml. ea.) 6 teal mixing tips	<u>Unflavored</u> #BR-9001 Regular Set #BR-9002 Fast Set #BR-9003 Super Fast Set	<u>Mint Chocolate Chip</u> #BR-9004 Fast Set #BR-9005 Super Fast Set <u>“HARD” Bite – Unflavored</u> #BR-9010 Fast Set
		BLU-BITE® HP Blu-Mousse®

*Legal disclaimers: BLU-BITE® is a registered trademark of Henry Schein®. Blu-Mousse® is a registered trademark of Parkell.

Use of these trademarks does not imply any affiliation with or endorsement by these companies. Mydent International is not responsible for any innacuracy or typographical errors in this document.

Request FREE Samples at: www.defend.com/freesamples

Mydent
International

80 Suffolk Court, Hauppauge, NY 11788
800-275-0020 • FAX 631-434-7750
www.defend.com • EMAIL: sales@defend.com



Follow us on Facebook.com/MydentDefend
Follow us on Twitter.com/Mydent_Defend