



**Works Better.
Lasts Longer.
Costs Less.**



EVACUATION



#DT-5501

#DT-5501

Fits A-dec models with black plastic canisters, plus newer models by Ampco, Belmont, DCI, Forest, Marus, MDT, & Schein - diameter 2 1/8".



#DT-5502

#DT-5502

Fits A-dec models with metal canisters - diameter 2 15/16".



#DT-5503

#DT-5503

Fits Ampco, Belmont, Biotec, Chayes, Dabi-Atlante, Dansereau, DCI, Den-Tal-EZ, Dentech, Forest, Knight, Pelton & Crane, Proma and others - diameter 1 7/8".

***Universal Diameter**



#DT-6200C

#DT-6200C*

Fits Siemens, Planmeca and many cuspidors with small diameter openings.



#DT-HVE

#DT-HVE*

Fits as a "cap" on all standard high volume evacuator tips.



#DT-5508S

#DT-5508S*

Fits all Saliva Ejectors.



#DT-5509

#DT-5509

Fits DCI/Belmont Reality models and System 6000 models - diameter 1 1/2"



#DT-6100C

#DT-6100C*

Fits most cuspidor bowls including: Adec, Ampco, Belmont, Dabi-Atlante, Dansereau, DCI, Den-Tal-EZ, Forest, Knight, Marus, MDT, Proma, Ritter & Schein.

"Works Better"

Sized to fit virtually every evacuation unit on the market today.

"Lasts Longer"

Provides the safety and convenience that you demand from a disposable infection control item.

"Costs Less"

The quality of a name brand item at the value of a private label.

Trap-EZ Disposable Traps

DEFEND Trap-EZ disposable filters are designed to replace the trap presently in your evacuation unit or cuspidor. They reduce the possibility of infection and contamination, while improving the performance of your evacuation system.

Packaging:	Item No:	Compares to*:
144 per Box	DT-5501 EVACUATION TRAP 2 1/8" DT-5502 EVACUATION TRAP 2 1/4" DT-5503 EVACUATION TRAP 1 7/8" DT-5508S SALIVA EJECTOR SCREEN DT-5509 EVACUATION TRAP 1 1/2" DT-6100C CUSPIDOR TRAP 144/BX DT-6200C CUSPIDOR TRAP UNIVERSAL WIDE FLANGE DT-HVE HVE TIP SCREEN	Dispos-A-Trap® - TotalCare featuring Pinnacle Products Dispos-A-Screen® Traps - Zirc

*Legal disclaimers: Dispos-A-Trap is a registered trademark of TotalCare. Dispos-A-Screen® traps is a registered trademark of Zirc. Use of these trademarks does not imply any affiliation with or endorsement by these companies. Mydent International is not responsible for any innaccuracy or typographical errors in this document.

Request FREE Samples at: www.defend.com/freesamples

Mydent
International

80 Suffolk Court, Hauppauge, NY 11788
800-275-0020 • FAX 631-434-7750
www.defend.com • EMAIL: sales@defend.com



Follow us on Facebook.com/MydentDefend
Follow us on Twitter.com/Mydent_Defend