




**Works Better.
Lasts Longer.
Costs Less.**



STERILIZATION

**“Works Better”
“Lasts Longer”
“Costs Less”**



Heat Sealer & Accessories

This ETL Listed, ultra durable heat sealer will seal Nylon Sterilization Tubing within 1-2 seconds and then turn off automatically. Strong stainless steel heating elements resists breaking.

Item No:	Description:	Compares to*:
HS-3000	110 Volt Unit	Nyclave® Impulse Heat sealer Impulse Heat Sealer - Henry Schein®
HS-3001	Heating Elements (3/pk)	
HS-3002PK	Upper Insulating Cloth (3/pk)	
HS-3003PK	Lower Insulating Cloth (3/pk)	
HS-3004	Silicone Strip (opaque plastic)	
HS-3005	Heat Sealer Repair Kit Contains: 1 Heating Element, 1 Upper Cloth & 1 Silicone Wrap	

*Legal disclaimers: Nyclave® is a registered trademark of Young Dental. Impulse Heat Sealer is a registered trademark of Henry Schein®. Use of these trademarks does not imply any affiliation with or endorsement by these companies. Mydent International is not responsible for any innacuracy or typographical errors in this document.



Nylon Sterilization Tubing

Manufactured from transparent nylon for easy viewing of contents, the sterilization tubing will remain pliable after heat sterilization. The tubing is for use in autoclave, dry heat or chemical sterilizers and may be sealed with a heat sealer. The tubing comes in a 100 foot roll and has an easy to read sterilization indicator.

Packaging:	Item No:	Compares to*:
100 foot roll with DEFEND Indicator Ink	AT-1018 -2" AT-1019- 3" AT-1020- 4" AT-1021- 6"	DH Nyclave® Dual-Purpose Nylon Sterilization Tubing - Henry Schein®

*Legal disclaimers: DH Nyclave® is a registered trademark of Young Dental. Dual-Purpose Nylon Sterilization Tubing is a registered trademark of Henry Schein®. Use of these trademarks does not imply any affiliation with or endorsement by these companies. Mydent International is not responsible for any innacuracy or typographical errors in this document.

Request FREE Samples at: www.defend.com/freesamples

**Mydent
International**

80 Suffolk Court, Hauppauge, NY 11788
800-275-0020 • FAX 631-434-7750
www.defend.com • EMAIL: sales@defend.com



Follow us on Facebook.com/MydentDefend
Follow us on Twitter.com/Mydent_Defend