



**Works Better.
Lasts Longer.
Costs Less.**



EVACUATION

**“Works Better”
“Lasts Longer”
“Costs Less”**

Standard Saliva Ejectors

Fit all standard saliva ejector hose ends. Pliable, holds their shape and the tip minimizes tissue impingement. Flexible "bonded" tips. Latex free.

Packaging:	Item No:	Compares to*:
100 per bag	#SE-7000 - white w/ white tip #SE-7001 - clear w/blue tip #SE-7002 - assorted w/ clear tip	Premium Saliva Ejectors - Crosstex Saliva Ejectors - Henry Schein®

*Legal disclaimers: Premium Saliva Ejectors is a registered trademark of Crosstex. Saliva Ejectors is a registered trademark of Henry Schein®.

Vented Oral Evacuation Tips

These evacuators fit all standard evacuator hose ends and contain smooth edges for maximum patient comfort. The straight design allows for increased control and they are convenient single-use disposable tips. One end is vented and the other is non-vented (dual use). Non sterile, 1/8" diameter.

Packaging:	Item No:	Compares to*:
100 per bag	#VE-1121 - 6" length	Combo Tip® Evacuator Tips Evacuator Tips - Henry Schein®

*Legal disclaimers: Combo Tip® is a registered trademark of Crosstex. Evacuator Tips is a registered trademark of Henry Schein®.

Surgical Aspirator Tips

Made with smooth edges for maximum patient comfort. The surgical aspirator is molded at a 30 degree angle and fits tall standard evacuator hose ends. This item is also autoclavable up to 135 C or may be disposed after single use.

Packaging:	Item No:	Compares to*:
25 per bag	#ST-1020 - blue 1/16" #ST-1021 - white 1/8" #ST-1023 - green 1/4"	Monoject™ Aspirator Tips #444 Surgical Aspirator Tips - Henry Schein®

*Legal disclaimers: Monoject™ Aspirator Tips #444 is a registered trademark of Covidien. Surgical Aspirator Tips is a registered trademark of Henry Schein®. Use of these trademarks does not imply any affiliation with or endorsement by these companies. Mydent International is not responsible for any innaccuracy or typographical errors in this document.

Request FREE Samples at: www.defend.com/freesamples

**Mydent
International**

80 Suffolk Court, Hauppauge, NY 11788
800-275-0020 • FAX 631-434-7750
www.defend.com • EMAIL: sales@defend.com



Follow us on Facebook.com/MydentDefend
Follow us on Twitter.com/Mydent_Defend