



**Works Better.  
Lasts Longer.  
Costs Less.**



## IMPRESSION



VPS Impression Materials Setting Times					
Viscosity		Color	Work Time	Set Time	Total Set Time
Light Body	Regular Set	Yellow	2:30	1:00	3:30
	Fast Set		1:00	1:00	2:00
Medium Body	Regular Set	Buff	2:30	1:00	3:30
	Fast Set		1:00	1:00	2:00
Monophase Body	Regular Set	Purple	2:30	1:00	3:30
	Fast Set		1:00	1:00	2:00
Heavy Body	Regular Set	Rose	2:30	1:00	3:30
	Fast Set		1:00	1:00	2:00

### “Works Better”

Contains wetting agents that enhance the hydrophilic properties during impression taking that allow the material to displace moisture, resulting in extremely accurate impressions.

### “Lasts Longer”

Excellent tear strength, will not deform or cause tissue inflammation allowing you to concentrate on taking the perfect impression.

### “Costs Less”

The superior formula saves on "re-takes" and is still about half the cost of the leading competitor.

## Super Hydrophilic Vinyl Polysiloxane Impression Material (VPS)

High performance (H.P) cartridge system available in regular and fast set in all 4 viscosities. The excellent tear strength allows the material to maintain "impression" integrity upon removal and precise impressions no matter what technique you use. Exceeds ISO 4380 standard which will help eliminate retakes.

Packaging:	Item No:	Compares to*:															
<b>4 cartridges</b> (50 ml. ea.)	<table> <tr> <td>Viscosity</td><td>Regular Set</td><td>Fast Set</td></tr> <tr> <td>Light</td><td>#VP-9003</td><td>#VP-9008</td></tr> <tr> <td>Medium</td><td>#VP-9004</td><td>#VP-9009</td></tr> <tr> <td>Mono</td><td>#VP-9005</td><td>#VP-9010</td></tr> <tr> <td>Heavy</td><td>#VP-9007</td><td>#VP-9011</td></tr> </table>	Viscosity	Regular Set	Fast Set	Light	#VP-9003	#VP-9008	Medium	#VP-9004	#VP-9009	Mono	#VP-9005	#VP-9010	Heavy	#VP-9007	#VP-9011	<b>VP MIX® High Performance Sharp™ Impression Material</b>
Viscosity	Regular Set	Fast Set															
Light	#VP-9003	#VP-9008															
Medium	#VP-9004	#VP-9009															
Mono	#VP-9005	#VP-9010															
Heavy	#VP-9007	#VP-9011															

\*Legal disclaimers: VP MIX is a registered trademark of Henry Schein®. Sharp™ Impression Material is a registered trademark of Parkell, Inc.

Use of these trademarks does not imply any affiliation with or endorsement by these companies. Mydent International is not responsible for any innaccuracy or typographical errors in this document.

**Request FREE Samples at: [www.defend.com/freesamples](http://www.defend.com/freesamples)**

**Mydent**  
International

80 Suffolk Court, Hauppauge, NY 11788  
800-275-0020 • FAX 631-434-7750  
[www.defend.com](http://www.defend.com) • EMAIL: [sales@defend.com](mailto:sales@defend.com)



Follow us on Facebook.com/MydentDefend  
Follow us on Twitter.com/Mydent\_Defend

