



**Works Better.
Lasts Longer.
Costs Less.**



DISPOSABLES



“Works Better”

Designed specifically for placing etchants, sealants, cements and other flowable materials.

“Lasts Longer”

The universal, color-coded, Luer-lock design makes it easy to match proper gauge to the corresponding material.

“Costs Less”

Delivers superior performance while providing comfort for the patient – all at an economical price.

Pre-Bent Dispensing Tips

These dispensing tips are pre-bent to the proper angle for placing etchants, sealants, cements and other flowable materials. German made & final step washed and blown dry to avoid sediment within the metal tubes. Packaged 100/bag and are color coded to represent gauges for each flowable material.

Packaging:	Item No:	Compares to*:
100 per bag	<p>#PN-1227 Pink [Liners & Bases (18 ga)]</p> <p>#PN-1228 Yellow [Sealant & Composites - 20ga.]</p> <p>#PN-1229 Black [Composites & Cements - 19ga.]</p> <p>#PN-1230 Grey [Sealant & Cement - 22ga.]</p> <p>#PN-1231 Blue [Etchants - 25ga.]</p>	Prebent Multipurpose Dispensing Tips - Henry Schein®

*Legal disclaimers: Prebent Multipurpose Dispensing Tips is a product of Henry Schein®.

Use of these trademarks does not imply any affiliation with or endorsement by these companies. Mydent International is not responsible for any innacuracy or typographical errors in this document.

Request FREE Samples at: www.defend.com/freesamples

Mydent
International

80 Suffolk Court, Hauppauge, NY 11788
800-275-0020 • FAX 631-434-7750
www.defend.com • EMAIL: sales@defend.com



Follow us on Facebook.com/MydentDefend
Follow us on Twitter.com/Mydent_Defend