


VACUUM FORMING EQUIPMENT & MATERIALS



Buffalo Dental Manufacturing Company, Inc.
www.buffalodental.com



The Dynamic Duo™

A time-saving knife with a hot-tipped blade that promotes fast and easy trimming of ALL vacuum formed materials. Adjustable temperature control for precise heat adjustment providing the operator with immediate heat for fast cutting and trimming. Ultra-fast heating circuit gets tip to ideal temperature in seconds, saving time. Lightweight pen with sharp, permanent tip promotes improved comfort and control. We have even added an EXTRA OUTLET for an additional Pen. Unit Comes complete with 1 Knife edged Tip Pen.

Item # Description

- 38760 Buffalo Dynamic Duo V-35/Therma-Knife Combo Unit Complete, 120 V AC
38765 Buffalo Dynamic Duo V-35/Therma-Knife Combo Unit Complete, 220 V AC



The Dynamic Duo™

V-35 MicroMotor – Electric Handpiece System

The V35 electric handpiece system features a no-frills, multi-speed control console teamed with our new Gray Handpiece. It represents a heavy-duty, maximum 35,000 RPM, lab-only handpiece system for use with HP-sized burs that offers high durability, high torque, and economical pricing.

Item # Description

- 37570 V35 Gray Electric Handpiece System, 120 V AC (V35 Console, Gray Electric Handpiece, Foot Control and Handpiece Rest)
37575 V35 Gray Electric Handpiece System, 220 V AC (V35 Console, Gray Electric Handpiece, Foot Control and Handpiece Rest)



V-35

ThermaKnife™ Hot Tipped Tray Trimming Knife

Effective and Time-Saving Knife For Fast and Easy Trimming of Bleaching Trays, Night Guards, Mouthguards, and Other Thermoplastic Materials.

- Fully adjustable temperature control for precise cutting action
- Cutting tip heats ultra-fast to save time and eliminate waiting
- Lightweight pen features an ergonomic design and sharp, permanent tip for improved comfort and control Knife-Edge Tip
- A & B selection switch allows for the use of two different Thermal Pens to reduce time and increase trimming productivity
- Fast and easy pen replacement (via snap connector at back of pen) for added convenience Concave Tip
- Convenient pen rest provides added safety when unit not in use
- Small portable size: Console 5.75" W x 5" D x 4"H (incl. Rest) Pen: 6" L x .75" W (at foam grip,) 2.5 lb. total console/pen weight

Item # Description

- 80500 ThermaKnife Intro Kit, 120V AC
(Includes Knife Edge Pen)

- 80510 Replacement Pen with Concave Tip
80520 Replacement Pen with Knife Edge Tip
80530 Replacement Pen with Waxing Tip



ThermaKnife w/ Knife Edge Tip

ThermaPlier™ Hot Tipped Thermal Pliers

For accentuating and making retentive undercuts on orthodontic appliances and mouthguards. Reheats the critical retention area allowing emphasis without causing adjacent plastic to distort.

- Adjust the temperature to a setting which allows the creation of dimples in the material you are working with but which does not burn the material.
- Use the plier to make retentive undercuts on the orthodontic appliances and mouthguards as per figure 1 & 2.




Figure 1



Figure 2



ThermaPlier™

Item # Description

- 80610 ThermaPlier 110V

Buffalo Dental Manufacturing Co., Inc.

159 Lafayette Drive • Syosset, NY 11791-0678 USA
Tel: (800) 828-0203 • (516) 496-7200 • Fax: (516) 496-7751 • Fax(orders): (516) 213-6610
www.buffalodental.com

Accu-Vac™ Precision Vacuum Forming System

An extremely powerful, accurate, and quiet vacuum forming system utilizing compressed air-powered suction.

- Utilizes ultra-high powered venturi suction (requires 60 to 100 psi external compressed air supply) to create highly accurate vacuum adaptations
- Gasket-free frame accepts 5" diameter round and 5" x 5" square materials for improved versatility

- Unique table uses either standard flat plate or included steel bead basket for multiple model positioning options
- Single On/Off heat control switch and automatic vacuum actuation switch promote easy operation
- Convenient small size (7 1/2" W x 10 1/2" D x 10" H x 16 lbs.) enhances portability
- System includes Accu-Vac unit, Flat Plate Table, Steel Beads, Bead Basket, and sample range of Vacuum Forming Materials

Item # Description

- | | |
|-------|---------------------------|
| 80175 | Accu-Vac System, 120 V AC |
| 80176 | Accu-Vac System, 220 V AC |



Accu-Vac

Sta-VacII™ Versatile Vacuum Forming System

Our improved classic high performance electric vacuum former that now allows both traditional flat table and highly accurate bead basket adaptations.

- New ConverTable™ combination vacuum table/bead basket can be used in standard flat table mode or reversed and used in bead basket mode for improved versatility, vacuum flow, and accuracy
- Dual-stage electric vacuum motor provides powerful suction and highly accurate adaptations- no compressed air connections required!
- Spring loaded, self-leveling, gasket-sealed frame handles a variety of 5" x 5" sheet materials, and the square center post assures smooth operation and guarantees alignment of frame and table
- Solid construction (7" W x 16" D x 11" H x 32 lbs.) for maximum durability and long life
- Unique housing design prevents motor contamination, shields heater, and promotes ergonomic operation

- All systems include a sample range of vacuum forming materials

Item # Description

- | | |
|-------|--|
| 80193 | Sta-Vac II Vacuum Forming System, 120 V AC, w/ Two ConverTables, Steel Beads and Plastic Work Tray |
| 80194 | Sta-Vac II Vacuum Forming System, 220 V AC, w/ Two ConverTables, Steel Beads and Plastic Work Tray |

Also available is the Sta-Vac, with all internal components identical to the Sta-Vac II, but with only one ConverTable and no steel beads nor work tray, designed for use in original flat model table mode only.

Item # Description

- | | |
|-------|--------------------------|
| 80191 | Sta-Vac System, 220 V AC |
| 80190 | Sta-Vac System, 120 V AC |



Sta-Vac II



Sta-Vac

Econo-Vac™ Vacuum Former

Premium vacuum former features power and performance in a standard, upright design.

- Dual-stage vacuum motor provides powerful suction and accurate adaptations
- Spring loaded, self-leveling, gasket-sealed frame handles a variety of 5" x 5" sheet materials

- Square center post assures smooth operation and guarantees perfect alignment of table and frame
- Solid construction (7" W x 10 1/2" D x 15" H x 32 lbs.) for high durability and long life
- All systems include a sample range of vacuum forming materials

Item # Description

- | | |
|-------|----------------------------|
| 80186 | Econo-Vac System, 120 V AC |
| 80187 | Econo-Vac System, 220 V AC |



Econo-Vac

Tray-Vac™ Vacuum Former

A high performance, no-frills, economical, and reliable small vacuum former for the dental office.

- High-velocity air flow system promotes rapid, highly accurate adaptations with 5" x 5" sheet materials
- Unique quick-lock frame and gasket-free design promote convenient, reliable operation

- Lightweight, small machine (7 1/2" W x 9 1/2" D x 15" H x 17 lbs.) improves portability
- All systems include a sample range of vacuum forming materials

Item # Description

- | | |
|-------|--------------------|
| 80165 | Tray-Vac, 120 V AC |
| 80168 | Tray-Vac, 220 V AC |



Tray-Vac



Sta-Lube™ Spray

A convenient silicone lubricant spray that prevents materials from sticking to models, vacuum table and frame.

- Environmentally safe material assures easy clean-up of flasks and trays
- Ideal separating medium for dies when making copings

Item # Description
80015 Sta-Lube Silicone Spray, 16 oz. Can



Sta-Lube

Alcohol Torch

An economical torch that produces a needle-point flame with only a gentle squeeze of the plastic bottle.

- 5" high torch uses denatured alcohol fuel for cost savings and convenience
- Base of bottle weighted with marbles to prevent tipping and improve safety
- Ideal for applying a rolled edge to vacuum formed adaptations (flame edge and rub with thumb)
- Can also be used as an alcohol lamp

Item # Description
82440 Alcohol Torch
57560-1D Wicks for Alcohol Torch, Pkg. of 12



Alcohol Torch

No. 55 Plate Shears

Heavy duty scissor-style instruments for multiple cutting and trimming applications.

- No. 61 plate shears excellent for trimming vacuum forming materials

Item # Description
78250 No. 55 Plate Shears



#55

Arbor Band Mandrels

High quality HP (3/32" diameter) mandrels with adjustable -diameter rubber heads are designed to firmly hold arbor bands in place during use.

- Replacement rubber heads available for renewing older mandrels

Item # Description
02771 3/8" Medium Arbor Bands , pkg 100
02920 3/8" Complete Mandrel, HP Sized Shank



Burs, Miscellaneous

High quality burs for miscellaneous applications.

- Cutting and Finishing Bur utilizes hardened tool steel for fast and accurate cutting to 15,000 rpm. Bur is ideal for trimming custom trays, base plates, and tough thermoplastic sheet materials.

- Plastic & Stone Trimming Burs designed for fast and smooth trimming of stone and epoxy dies, models, acrylics, and sheet resin materials. Burs constructed of carbide steel for long life. Popular taper shape with non-clogging flutes.

Item # Description
20000 Cutting & Finishing Bur for Thermoplastic Materials, HP



Cutting & Finishing Bur

PRKIT

Item # Description
18100 Cone 82T carbide bur
18010 Taper 84T carbide bur
02781 3/8" Coarse arbor bands (100)
02920 3/8" Complete arbor band mandrel

Block-Out Compound

Used for blocking out remaining teeth and undercuts on models when making appliances.

- Reusable, heat resistant material resists running and sticking

- Can also be rolled around the base of models to promote easy periphery trimming

Item # Description
05240 Block-Out Material, 1 lb.

Buffalo Dental Manufacturing Co., Inc.

159 Lafayette Drive • Syosset, NY 11791-0678 USA

Tel: (800) 828-0203 • (516) 496-7200 • Fax: (516) 496-7751 • Fax(orders): (516) 213-6610

www.buffalodental.com



Vacuum Forming Materials

High quality 5" x 5" sheet thermoplastic materials for creating a wide range of vacuum formed appliances and trays.

- Pure grade resin sheets are latex-free and uniform in size, cut and thickness for improved performance.
- Resist bubbling, crazing and fracturing to promote superior adaptations every time.
- Clear soft Ethyl Vinyl Acetate (EVA) sheets designed to bond together for build-up and lamination of custom mouthguards.

Rigid Clear Tray and Splint Material

Buffalo Clear Material offers the best model adaptation and excellent rigidity. The result will be a clear, highly polished surface with no color dyes or bluing agents. The Extra Wide size stays secure in all frames.

.020" (0.5mm) Clear Tray Material

| Item # | Description |
|--------|--|
| 81660 | Pkg. of 50 81665 Pkg. of 200 81670 Pkg. of 625 |

.040" (1.0mm) Clear Tray Material

| Item # | Description |
|--------|-----------------------------------|
| 30990 | Pkg. of 25 30995 Pkg. of 625 |

.060" (1.5mm) Clear Tray Material

| Item # | Description |
|--------|-----------------------------------|
| 31000 | Pkg. of 25 31010 Pkg. of 625 |

.080" (2.0mm) Clear Tray Material

| Item # | Description |
|--------|-----------------------------------|
| 31090 | Pkg. of 25 31100 Pkg. of 625 |

.100" (2.5mm) Clear Tray Material

| Item # | Description |
|--------|-----------------------------------|
| 31110 | Pkg. of 25 31120 Pkg. of 625 |

.020" Coping and Temporary C&B Material

Cloudy, stiff Polypropylene sheets excellent for use when fabricating copings. Sized for improved vacuum is easier to use and promotes improved vacuum seal. Will not adhere to acrylics.

| Item # | Description |
|--------|--|
| 31720 | Pkg. of 50 31725 Pkg. of 200 31730 Pkg. of 625 |

.060" Pink Base Plate Material

Stiff pink Polystyrene sheets excellent for fabrication of thin base plates and custom trays.

| Item # | Description |
|--------|--|
| 04610 | Pkg. of 25 04615 Pkg. of 200 04620 Pkg. of 625 |

.080" Pink Tray Material

Stiff pink Polystyrene sheets excellent for fabrication of medium and thick base plates and custom trays.

| Item # | Description |
|--------|--|
| 82610 | Pkg. of 25 82615 Pkg. of 200 82620 Pkg. of 625 |

.100" Pink Tray Material

Stiff pink Polystyrene sheets excellent for fabrication of thick base plates and custom trays.

| Item # | Description |
|--------|--|
| 82650 | Pkg. of 25 82655 Pkg. of 200 82660 Pkg. of 625 |

.125" Blue Tray Material

Stiff blue Polystyrene sheets ideal for fabrication of custom trays.

| Item # | Description |
|--------|--|
| 82690 | Pkg. of 25 82695 Pkg. of 200 82700 Pkg. of 625 |

.150" Blue Tray Material

Stiff blue Polystyrene sheets ideal for fabrication of thick custom trays.

| Item # | Description |
|--------|--|
| 82730 | Pkg. of 25 82735 Pkg. of 200 82740 Pkg. of 625 |

.040" Clear Bleaching Tray Material

Clear, soft EVA Sheets excellent for use when fabricating bleaching trays. Sized for improved vacuum is easier to use and promotes improved vacuum seal.

| Item # | Description |
|--------|--|
| 62820 | Pkg. of 12 62803 Pkg. of 25 62830 Pkg. of 100 62840 Pkg. of 300 |

.080" Clear Thin Mouthguard Material

Clear, soft EVA Sheets excellent for fabrication of thin athletic mouthguards for non-contact sports, bruxism splints, soft splints, and custom fluoride trays. Material bonds to itself so it can be built up in layers for creation of athletic mouthguards.

| Item # | Description |
|--------|--|
| 62804 | Pkg. of 12 62806 Pkg. of 25 62802 Pkg. of 100 62800 Pkg. of 300 |

.118" Clear Mouthguard Material

Clear, soft EVA Sheets excellent for fabrication of light athletic mouthguards for non-contact sports, soft splints, and custom fluoride trays. Material bonds to itself so it can be built up in layers for creation of athletic mouthguards.

| Item # | Description |
|--------|--|
| 62815 | Pkg. of 12 62816 Pkg. of 25 62805 Pkg. of 100 62810 Pkg. of 300 |

.160" Clear Mouthguard Material

Thick, clear, resilient EVA sheets for fabrication of thick athletic mouthguards and resilient bruxism splints. Material bonds to thin EVA mouthguard materials so adaptations can be built up in layers for creation of athletic mouthguards.

| Item # | Description |
|--------|--|
| 62700 | Pkg. of 12 62725 Pkg. of 25 62705 Pkg. of 100 62710 Pkg. of 300 |

.160" Sta-Guard™ Clear Mouthguard Material

Thick, clear, resilient EVA classic-brand sheets excellent for fabrication of thick athletic mouthguards and resilient bruxism splints. Material bonds to thin EVA mouthguard materials so adaptations can be built up in layers for creation of athletic mouthguards.

| Item # | Description |
|--------|--|
| 62750 | Pkg. of 12 62755 Pkg. of 100 62760 Pkg. of 300 |



Thick resilient 5" X 5" EVA Colored Mouthguard Sheets, Excellent Protection for All Sports!

Buffalo Dental Mouthguard & Tray material feature a "Total Smooth" surface for comfort and appearance! Contains *NO BPA, NO PVC, NO LATEX, or PVC INSERTS!* Also features NO inclusions or glitter which can be swallowed! This material is used for athletic mouthguards.

.040" EVA Mouthguard Material 62610 PINK, Pkg. of 12

Item # Description

| | |
|-------|---------------------|
| 72560 | RED, Pkg. of 12 |
| 72570 | RED, Pkg. of 100 |
| 72580 | RED, Pkg. of 300 |
| 72590 | BLUE, Pkg. of 12 |
| 72592 | BLUE, Pkg. of 100 |
| 72593 | BLUE, Pkg. of 300 |
| 72595 | ORANGE, Pkg. of 12 |
| 72596 | ORANGE, Pkg. of 100 |
| 72597 | ORANGE, Pkg. of 300 |
| 72585 | BLACK, Pkg. of 12 |
| 72586 | BLACK, Pkg. of 100 |
| 72587 | BLACK, Pkg. of 300 |
| 72581 | WHITE, Pkg. of 12 |
| 72582 | WHITE, Pkg. of 100 |
| 72583 | WHITE, Pkg. of 300 |
| 72476 | YELLOW, Pkg. of 12 |
| 72477 | YELLOW, Pkg. of 100 |
| 72478 | YELLOW, Pkg. of 300 |

.040" EVA Mouthguard Material Disc (120mm)

Item # Description

| | |
|-------|-------------------|
| 62827 | CLEAR, Pkg. of 12 |
|-------|-------------------|

.040" EVA Mouthguard Material Disc (125mm)

Item # Description

| | |
|-------|--------------------|
| 62827 | CLEAR, Pkg. of 12 |
| 62822 | CLEAR, Pkg. of 100 |
| 62825 | CLEAR, Pkg. of 300 |
| 72670 | RED, Pkg. of 12 |
| 72680 | RED, Pkg. of 100 |
| 72690 | RED, Pkg. of 300 |
| 72691 | BLACK, Pkg. of 12 |
| 72692 | BLACK, Pkg. of 100 |
| 72693 | BLACK, Pkg. of 300 |



MADE IN USA

.080" EVA Mouthguard Material

Item # Description

| | |
|-------|-------------------------|
| 62788 | BLACK, Pkg. of 12 |
| 62795 | BLACK, Pkg. of 300 |
| 62789 | WHITE, Pkg. of 12 |
| 62799 | WHITE, Pkg. of 300 |
| 62784 | ORANGE, Pkg. of 12 |
| 62798 | ORANGE, Pkg. of 300 |
| 62792 | YELLOW, Pkg. of 12 |
| 62794 | YELLOW, Pkg. of 300 |
| 62797 | RED, Pkg. of 12 |
| 62796 | RED, Pkg. of 300 |
| 62670 | ROYAL BLUE, Pkg. of 12 |
| 62680 | ROYAL BLUE, Pkg. of 100 |
| 62690 | ROYAL BLUE, Pkg. of 300 |
| 62655 | PINK, Pkg. of 12 |
| 62664 | PINK, Pkg. of 100 |
| 62675 | PINK, Pkg. of 300 |

.080" EVA Mouthguard Material Disc (120mm)

Item # Description

| | |
|-------|-------------------------|
| 62814 | CLEAR, Pkg. of 100 |
| 62812 | CLEAR, Pkg. of 300 |
| 62671 | ROYAL BLUE, Pkg. of 12 |
| 62681 | ROYAL BLUE, Pkg. of 100 |
| 62691 | ROYAL BLUE, Pkg. of 300 |
| 62666 | PINK, Pkg. of 12 |
| 62676 | PINK, Pkg. of 100 |
| 62686 | PINK, Pkg. of 300 |
| 62791 | RED, Pkg. of 300 |
| 62793 | YELLOW, Pkg. of 300 |
| 62823 | BLACK, Pkg. of 300 |

.080" EVA Mouthguard Material Disc (125mm)

Item # Description

| | |
|-------|---------------------|
| 62807 | CLEAR, Pkg. of 12 |
| 62808 | CLEAR, Pkg. of 100 |
| 62809 | CLEAR, Pkg. of 300 |
| 62813 | WHITE, Pkg. of 300 |
| 62824 | BLACK, Pkg. of 300 |
| 62835 | ORANGE, Pkg. of 300 |

.118" EVA Mouthguard Material

Item # Description

| | |
|-------|-------------------------|
| 72530 | WHITE, Pkg. of 12 |
| 72540 | WHITE, Pkg. of 100 |
| 72550 | WHITE, Pkg. of 300 |
| 72532 | YELLOW, Pkg. of 12 |
| 72534 | BLUE, Pkg. of 12 |
| 72536 | ORANGE, Pkg. of 12 |
| 72539 | ORANGE, Pkg. of 300 |
| 72538 | RED, Pkg. of 12 |
| 72537 | RED, Pkg. of 300 |
| 72500 | BLACK, Pkg. of 12 |
| 72510 | BLACK, Pkg. of 100 |
| 72520 | BLACK, Pkg. of 300 |
| 62640 | ROYAL BLUE, Pkg. of 12 |
| 62641 | ROYAL BLUE, Pkg. of 100 |
| 62642 | ROYAL BLUE, Pkg. of 300 |

.118" EVA Mouthguard Material Discs (120mm)

Item # Description

| | |
|-------|-------------------------|
| 62889 | BLUE, Pkg. of 300 |
| 62888 | CLEAR, Pkg. of 300 |
| 62890 | RED, Pkg. of 300 |
| 62891 | BLACK, Pkg. of 300 |
| 62899 | WHITE, Pkg. of 300 |
| 62895 | YELLOW, Pkg. of 100 |
| 62898 | YELLOW, Pkg. of 300 |
| 62643 | ROYAL BLUE, Pkg. of 300 |
| 62644 | ROYAL BLUE, Pkg. of 300 |
| 62645 | ROYAL BLUE, Pkg. of 300 |

.118" EVA Mouthguard Material Discs (125mm)

Item # Description

| | |
|-------|-----------------------|
| 62818 | CLEAR, Pkg. of 12 |
| 62817 | CLEAR, Pkg. of 100 |
| 62819 | CLEAR, Pkg. of 300 |
| 62887 | SKY-BLUE, Pkg. of 300 |
| 72630 | BLACK, Pkg. of 12 |
| 72640 | BLACK, Pkg. of 100 |
| 72650 | BLACK, Pkg. of 300 |
| 72600 | WHITE, Pkg. of 12 |
| 72610 | WHITE, Pkg. of 100 |
| 72660 | WHITE, Pkg. of 300 |
| 72620 | BLUE, Pkg. of 300 |

Buffalo Dental Manufacturing Co., Inc.

159 Lafayette Drive • Syosset, NY 11791-0678 USA

Tel: (800) 828-0203 • (516) 496-7200 • Fax: (516) 496-7751 • Fax(orders): (516) 213-6610

www.buffalodental.com



EVA Colored Mouthguard Sheets continued...

.160" EVA Mouthguard Material

| Item # | Description |
|--------|-------------------------|
| 62770 | YELLOW, Pkg. of 12 |
| 62775 | YELLOW, Pkg. of 25 |
| 62771 | BLUE, Pkg. of 12 |
| 62777 | BLUE, Pkg. of 25 |
| 62720 | SKY-BLUE, Pkg. of 12 |
| 62740 | SKY-BLUE, Pkg. of 100 |
| 62730 | SKY-BLUE, Pkg. of 300 |
| 62772 | RED, Pkg. of 12 |
| 62778 | RED, Pkg. of 25 |
| 62773 | ORANGE, Pkg. of 12 |
| 62779 | ORANGE, Pkg. of 25 |
| 62774 | BLACK, Pkg. of 12 |
| 62781 | BLACK, Pkg. of 25 |
| 62783 | BLACK, Pkg. of 300 |
| 62776 | WHITE, Pkg. of 12 |
| 62782 | WHITE, Pkg. of 25 |
| 62700 | CLEAR, Pkg. of 12 |
| 62725 | CLEAR, Pkg. of 25 |
| 62600 | PINK, Pkg. of 12 |
| 62630 | PINK, Pkg. of 100 |
| 62650 | PINK, Pkg. of 300 |
| 626100 | ROYAL BLUE, Pkg. of 12 |
| 626105 | ROYAL BLUE, Pkg. of 100 |
| 626110 | ROYAL BLUE, Pkg. of 300 |

.160" EVA Mouthguard Material Disc (125mm)

| Item # | Description |
|--------|--------------------|
| 62715 | CLEAR, Pkg. of 12 |
| 62718 | CLEAR, Pkg. of 100 |
| 62719 | CLEAR, Pkg. of 300 |



MADE IN USA

BUFF-TUFF Mouthguard Laminate

Buff-Tuff Mouthguard Laminate features cross-graining for extra strength and constant thickness of finished mouthguards. Injection molded mouthguards cannot be cross grained! Contains *NO BPA, DIOXIN, PVC or INTERNAL BRACES!* Will not separate! **Designed Specifically For Athletic use.**

BUFF-TUFF Mouthguard Laminate Squares & Discs

| Item # | Description |
|--------|---|
| 70000 | BUFF-TUFF Clear 5"x5", Square, Pkg. of 10 |
| 70010 | BUFF-TUFF Clear 5"x5", Square, Pkg. of 100 |
| 70100 | BUFF-TUFF Clear 5"x5", Disc, Pkg. of 10 |
| 70110 | BUFF-TUFF Clear 5"x5", Disc, Pkg. of 100 |
| 70120 | BUFF-TUFF Blue 5"x5", Square, Pkg. of 10 |
| 70130 | BUFF-TUFF Blue 5"x5", Square, Pkg. of 100 |
| 70140 | BUFF-TUFF Blue 5"x5", Disc, Pkg. of 10 |
| 70150 | BUFF-TUFF Blue 5"x5", Disc, Pkg. of 100 |
| 70160 | BUFF-TUFF Orange 5"x5", Square, Pkg. of 10 |
| 70170 | BUFF-TUFF Orange 5"x5", Square, Pkg. of 100 |
| 70180 | BUFF-TUFF Yellow 5"x5", Square, Pkg. of 10 |
| 70190 | BUFF-TUFF Yellow 5"x5", Square, Pkg. of 100 |
| 70200 | BUFF-TUFF Black 5"x5", Square, Pkg. of 10 |



| Item # | Description |
|--------|---|
| 70210 | BUFF-TUFF Black 5"x5", Square, Pkg. of 100 |
| 70220 | BUFF-TUFF Red 5"x5", Square, Pkg. of 10 |
| 70230 | BUFF-TUFF Red 5"x5", Square, Pkg. of 100 |
| 70240 | BUFF-TUFF Sky-Blue 5"x5", Square, Pkg. of 10 |
| 70250 | BUFF-TUFF Sky-Blue 5"x5", Square, Pkg. of 100 |
| 70260 | BUFF-TUFF White 5"x5", Square, Pkg. of 10 |
| 70270 | BUFF-TUFF White 5"x5", Square, Pkg. of 100 |

BUFF-TUFF Mouthguard Laminate Discs (125mm)

| Item # | Description | Item # | Description |
|--------|---|--------|---|
| 70280 | BUFF-TUFF Orange 5"x5", Disc, Pkg. of 10 | 70350 | BUFF-TUFF Red 5"x5", Disc, Pkg. of 100 |
| 70290 | BUFF-TUFF Orange 5"x5", Disc, Pkg. of 100 | 70360 | BUFF-TUFF Sky-Blue 5"x5", Disc, Pkg. of 10 |
| 70300 | BUFF-TUFF Yellow 5"x5", Disc, Pkg. of 10 | 70370 | BUFF-TUFF Sky-Blue 5"x5", Disc, Pkg. of 100 |
| 70310 | BUFF-TUFF Yellow 5"x5", Disc, Pkg. of 100 | 70380 | BUFF-TUFF White 5"x5", Disc, Pkg. of 10 |
| 70320 | BUFF-TUFF Black 5"x5", Disc, Pkg. of 10 | 70390 | BUFF-TUFF White 5"x5", Disc, Pkg. of 100 |
| 70330 | BUFF-TUFF Black 5"x5", Disc, Pkg. of 100 | 70400 | BUFF-TUFF Clear 5"x5", Disc, Pkg. of 10 |
| 70340 | BUFF-TUFF Red 5"x5", Disc, Pkg. of 10 | 70410 | BUFF-TUFF Clear 5"x5", Disc, Pkg. of 100 |

Buffalo Dental Manufacturing Co., Inc.

159 Lafayette Drive • Syosset, NY 11791-0678 USA

Tel: (800) 828-0203 • (516) 496-7200 • Fax: (516) 496-7751 • Fax(orders): (516) 213-6610

www.buffalodental.com



Colored Mouthguard Materials with Polyshok™ Technology!

New Polyshok™ technology allows the creation of high-impact sports mouthguards that can absorb 150% more energy than standard EVA mouthguards. The improved shock absorption can result in less destructive forces being transferred to dental structures, fewer oral-fascial injuries, and less concussions.

Item # Description

62901 .125" Polyshok™ Variety Pack, assorted box of 38 units

Item # Description

62931 .155" Polyshok™ Variety Pack, assorted box of 30 units

Item # Description

62902 .125" Black- Pkg. of 18

62903 .125" Blue- Pkg. of 18

62904 .125" Red- Pkg. of 18

62905 .125" Yellow- Pkg. of 18

62906 .125" Orange- Pkg. of 18

62907 .125" Purple- Pkg. of 18

62908 .125" Green- Pkg. of 18

62909 .125" White- Pkg. of 18

Item # Description

62902 .125" Black- Pkg. of 18

62903 .125" Blue- Pkg. of 18

62904 .125" Red- Pkg. of 18

62905 .125" Yellow- Pkg. of 18

62906 .125" Orange- Pkg. of 18

62907 .125" Purple- Pkg. of 18

62908 .125" Green- Pkg. of 18

62909 .125" White- Pkg. of 18



Polyshok™ GAMEguard (Pre-Formed "Boil & Bite")

High Performance Mouthguards w/ Carrying Cases- New! Polyshok™ GAMEguard SSD (No Strap) w/Case- Sports-Specific Design- Optimized Protection for Non-Helmeted Sports.

Item # Description

62970 Large Red SSD Pre-Formed Mouthguard

62971 Large Blue SSD Pre-Formed Mouthguard

62972 Large Black SSD Pre-Formed Mouthguard

62973 Large Neon Green SSD Pre-Formed Mouthguard

62974 Large Neon Yellow SSD Pre-Formed Mouthguard

62975 Small Red SSD Pre-Formed Mouthguard

62976 Small Blue SSD Pre-Formed Mouthguard

62977 Small Black SSD Pre-Formed Mouthguard

62978 Small Neon Green SSD Pre-Formed Mouthguard

62979 Small Neon Yellow SSD Pre-Formed Mouthguard



Large Red SSD Pre-Formed Mouthguard

Polyshok™ GAMEguard PLUS w/ Strap & Case- Optimized Extra Protection for Helmeted Sports.

Item # Description

62980 Large Red PLUS Pre-Formed Mouthguard

62981 Large Blue PLUS Pre-Formed Mouthguard

62982 Large Black PLUS Pre-Formed Mouthguard

62983 Large Green PLUS Pre-Formed Mouthguard

62984 Large Yellow PLUS Pre-Formed Mouthguard

62985 Small Red PLUS Pre-Formed Mouthguard

62986 Small Blue PLUS Pre-Formed Mouthguard

62987 Small Black PLUS Pre-Formed Mouthguard

62988 Small Green PLUS Pre-Formed Mouthguard

62989 Small Yellow PLUS Pre-Formed Mouthguard



Large Red PLUS Pre-Formed Mouthguard

Buffalo Dental Manufacturing Co., Inc.

159 Lafayette Drive • Syosset, NY 11791-0678 USA

Tel: (800) 828-0203 • (516) 496-7200 • Fax: (516) 496-7751 • Fax(orders): (516) 213-6610

www.buffalodental.com