




**Works Better.
Lasts Longer.
Costs Less.**



DISPOSABLES



Superior Quality & Less Expensive Replacement for:

Sure Tip – Henry Schein
Crystal Tip – Westside Resources
FASTips – Dentsply Professional
Safe-Tips EZ – Totalcare / Pinnacle
Spree – Totalcare / Pinnacle

“Works Better”

The durable, all plastic construction provides dry air every time with no water leakage.

“Lasts Longer”

No adaptors or conversion kits needed and will fit most models of equipment.

“Costs Less”

Delivers superior air/water spray at a price that cannot be beat.

Air/Water 3-Way Syringe Tips

DEFEND Air/Water 3-Way Syringe Tips are an excellent alternative to typical metal autoclavable tips. The disposable tips have an inner core water chamber that seats into the syringe, eliminating water contamination in the air chamber. The unique quad core design ensures instant dry air, prevents cross contamination, and delivers superior air/water spray. DEFEND syringe tips feature an O-ring groove that secures tip lock for patient safety and eliminates leakage. The patient angle is identical to that of metal tips, providing ideal placement of the tip anywhere in a patient’s mouth, and is strong enough to allow for retraction. More cost effective than comparable disposable tips.

Packaging:	Item No:	Compares to*:
250 per bag	#AW-1000 Clear	Sure Tip®
	#AW-2000 White	Flash Tips™
	#AW-3000 Blue	
1500 per bag	#AW-3500 Blue (Bulk Pack)	

*Legal disclaimers: Sure Tip® Plus is a registered trademark of Henry Schein®. Flash Tips™ is a registered trademark of Sultan Healthcare. Use of these trademarks does not imply any affiliation with or endorsement by these companies. Mydent International is not responsible for any inaccuracy or typographical errors in this document.

Request FREE Samples at: www.defend.com/freesamples

**Mydent
International**

80 Suffolk Court, Hauppauge, NY 11788
800-275-0020 • FAX 631-434-7750
www.defend.com • EMAIL: sales@defend.com



Follow us on Facebook.com/MydentDefend
Follow us on Twitter.com/Mydent_Defend

