




**Works Better.  
Lasts Longer.  
Costs Less.**



## DISPOSABLES



### “Works Better”

Siliconized to help ensure a gentle insertion of needle into tissue for minimal discomfort.

### “Lasts Longer”

The Prethreaded plastic hub for simple, straight attachment to the syringe perfectly and securely every time.

### “Costs Less”

Delivers superior performance while providing comfort for the patient – all at an economical price.

## Plastic Hub Self-Threading Dental Needles

Needle made of stainless steel; polypropylene hub with aluminum insert.

- Lot #, needle size and sterile seal on each individual needle casing.
- Disposable, sterile needles for use on standard 1.8mL dental syringes
- Triple-beveled needle features an arrow indicator to help orient the bevel position during injections
- Siliconized to help ensure a gentle insertion of needle into tissue for minimal discomfort
- Prethreaded plastic hub for simple, straight attachment to the syringe perfectly and securely everytime
- Color coding for easy needle diameter identification

Packaging:	Item No:	Compares to*:
100 per box	<b>#DN-2000 – 27GA LONG</b> <b>#DN-2100 – 27GA SHORT</b> <b>#DN-3000 – 30GA SHORT</b> <b>#DN-3100 – 30GA EXTRA SHORT</b>	<b>Premium Dental Needles - Plastic Hub - HenrySchein®</b>

\*Legal disclaimers: Premium Dental Needles - Plastic Hubs is a product of Henry Schein®.

Use of these trademarks does not imply any affiliation with or endorsement by these companies. Mydent International is not responsible for any innacuracy or typographical errors in this document.

**Request FREE Samples at: [www.defend.com/freesamples](http://www.defend.com/freesamples)**

**Mydent**  
International

80 Suffolk Court, Hauppauge, NY 11788  
800-275-0020 • FAX 631-434-7750  
[www.defend.com](http://www.defend.com) • EMAIL: [sales@defend.com](mailto:sales@defend.com)



Follow us on Facebook.com/MydentDefend  
Follow us on Twitter.com/Mydent\_Defend