



**Works Better.  
Lasts Longer.  
Costs Less.**



## IMPRESSION



### “Works Better”

The same benefits as our cartridge delivery VPS material in a hand mix option.

### “Lasts Longer”

Technique sensitive to provide the results that you depend upon.

### “Costs Less”

The most economical way to use VPS material.

VPS Impression Materials Setting Times					
Item#	Viscosity	Color	Work Time	Set Time	Total Set Time
VP-8006	Putty	Regular Set	2:00	2:00	4:00
VP-8008		Fast Set	1:00	2:00	3:00

## Super Hydrophilic Vinyl Polysiloxane Putty Material (VPS)

Hand Mix Extra Heavy Body VPS available in fast and regular set. The excellent tear strength allows the material to maintain "impression" integrity upon removal and precise impressions no matter what technique you use. Exceeds ISO 4380 standard which will help eliminate retakes.

Packaging:	Item No:	Compares to*:
2x jars 300 ml each of base & catalyst (600ml total)	VP-8006- Regular Set VP-8008- Fast Set	VP MIX®

\*Legal disclaimers: VP MIX is a registered trademark of Henry Schein®. Use of these trademarks does not imply any affiliation with or endorsement by these companies. Mydent International is not responsible for any innaccuracy or typographical errors in this document.

**Request FREE Samples at: [www.defend.com/freesamples](http://www.defend.com/freesamples)**

**Mydent**  
International

80 Suffolk Court, Hauppauge, NY 11788  
800-275-0020 • FAX 631-434-7750  
[www.defend.com](http://www.defend.com) • EMAIL: [sales@defend.com](mailto:sales@defend.com)



Follow us on Facebook.com/MydentDefend  
Follow us on Twitter.com/Mydent\_Defend

