



**Works Better.  
Lasts Longer.  
Costs Less.**



## IMPRESSION



### “Works Better”

Creates excellent temporary crowns, bridges, inlays/onlays and veneers with ease.

### “Lasts Longer”

Low shrinkage, low heat and a high compressive strength combine in a fast-setting provisional material that is second to none.

### “Costs Less”

The easy dispensing, MMA-free material provides consistent results that stack up to even the most expensive products on the market today.

## Temporary Crown and Bridge Material

DEFEND Temporary Crown & Bridge Material is a bis-acryl composite for chairside provisional restorations. DEFEND Temporary Crown & Bridge is fast and accurate, especially when used in conjunction with a quality vinyl polysiloxane impression material such as DEFEND Super Hydrophilic VPS (Monophase). Available in 6 shades: A1, A2, A3, A3.5, B1, and Bleach BL. All are delivered in 76g (50mL) cartridges designed to fit on a 10:1 style automix gun.

Packaging:	Item No:	Compares to*:
<b>1 76gm 10:1 Cartridge 10 mixing tips</b>	<b>CB-9000 A1      CB-9003 B1</b> <b>CB-9001 A2      CB-9004 BLEACH</b> <b>CB-9002 A3.5    CB-9005 A3</b>	<b>MAXITEMP® Luxatemp®</b>

\*Legal disclaimers: MAXITEMP® is a registered trademark of Henry Schein®. LUXATEMP® is a registered trademark of DMG America. Use of these trademarks does not imply any affiliation with or endorsement by these companies. Mydent International is not responsible for any innaccuracy or typographical errors in this document.

**Mydent  
International**

80 Suffolk Court, Hauppauge, NY 11788  
800-275-0020 • FAX 631-434-7750  
www.defend.com • EMAIL: sales@defend.com



Follow us on Facebook.com/MydentDefend  
Follow us on Twitter.com/Mydent\_Defend