



**Works Better.
Lasts Longer.
Costs Less.**



IMPRESSION MATERIAL ACCESSORIES

**“Works Better”
“Lasts Longer”
“Costs Less”**



Vinyl Tray Adhesive

This fast-acting tray adhesive provides firm adhesion between impression trays and VPS impression materials. This adhesive is compatible with virtually every “VPS” material on the market. Available in a 10ML w/ applicator.

Packaging:	Item No:	Compares to*:
10ml	VP-8300	VPS Tray Adhesive

*Legal disclaimers: VPS Tray Adhesive is a registered trademark of Henry Schein®. Use of these trademarks does not imply any affiliation with or endorsement by these companies. Mydent International is not responsible for any innacuracy or typographical errors in this document.



Standard Dispensing Gun

The dispensing gun is used to extrude impression materials, bite registration materials and temporary crown & bridge materials from HP cartridges. The durable release lever and plunger, which are coordinated simultaneously to load and remove the cartridge work consistently, with the trigger, to dispense the material smoothly and precisely. Works with 1 to 1 and 2 to 1 50ml cartridges.

Packaging:	Item No:	Compares to*:
1each	VP-8202	Dispensing Gun & Extruder™

*Legal disclaimers: Dispensing Gun is a registered trademark of Henry Schein®. Extruder™ is a registered trademark of Kerr™. Use of these trademarks does not imply any affiliation with or endorsement by these companies. Mydent International is not responsible for any innacuracy or typographical errors in this document.



4:1/10:1 Dispensing Gun

The dispensing gun is used to extrude impression materials and bite registration materials from cartridges. It consists of a release lever and a plunger, which are coordinated simultaneously to load and remove the cartridge. The gun also consists of a trigger to dispense the material. Works with 4 to 1 and 10 to 1 50ml cartridges.

Packaging:	Item No:	Compares to*:
1each	VP-8204	Dispensing Gun & Extruder™

*Legal disclaimers: Dispensing Gun is a registered trademark of Henry Schein®. Extruder™ is a registered trademark of Kerr™. Use of these trademarks does not imply any affiliation with or endorsement by these companies. Mydent International is not responsible for any innacuracy or typographical errors in this document.

Request FREE Samples at: www.defend.com/freesamples

**Mydent
International**

80 Suffolk Court, Hauppauge, NY 11788
800-275-0020 • FAX 631-434-7750
www.defend.com • EMAIL: sales@defend.com



Follow us on Facebook.com/MydentDefend
Follow us on Twitter.com/Mydent_Defend