



ProEdge
DENTAL WATER LABS

THEIR PROTECTION
YOUR ASSURANCE

*We leverage the full expertise
of our labs to protect your
patients, practice, and reputation.*



ProEdgeDental.com

Toll Free - 888-THE-EDGE (843.3343)

7042 South Revere Parkway, Suite 400
Centennial, Colorado 80112

DENTAL WATERLINE TESTING

*Because
you care
about them
and we care
about you*

The only way to know your practice meets the CDC and ADA guidelines is to regularly test your water.

Too many offices risk their practice and their reputation by not testing. Another 1/3 of those who do a test, still fail to meet the 500 CFU/ml guideline. Why? Because waterline maintenance is not easy and some products don't work very well.

ProEdge Dental Water Labs has professionally tested dental waterlines well over 100,000 times. 14 years experience testing dental water and providing sound, scientifically based advice has taught us how to best help practices around the country find a safe water solution.

We like the science, but what we care about is you getting back to what you do best - caring for patients.

72


children

were hospitalized in 2016 with serious bacterial infections after routine baby root canals (pulpotomies). The health department confirmed that contaminated dental waterlines were the cause.



HOW DOES YOUR PRODUCT PERFORM?

Pass Rate for Waterline Treatments



*See specific product results at ProEdgeDental.com

THE CDC'S STANDARD OF CARE

All Dental settings should:

1. Use water that meets the 500CFU/mL standard for routine procedures
2. Treat Dental Unit Water Lines (using distilled or untreated water is insufficient)
3. Monitor frequently per manufacturers recommendations
4. Use sterile water or saline for surgery

*Source: <https://www.cdc.gov/oralhealth/infectioncontrol/pdf/safe-care2.pdf>

TEST WITH PROEDGE

A Simple Mail-In Testing Process:



1. TAKE A WATER
SAMPLE DIRECT
FROM WATERLINE



2. SHIP OVERNIGHT
TO PROEDGE LABS



3. PROEDGE LABS
PROCESSES
THE SAMPLE



4. WE EMAIL YOUR
RESULTS REPORT



5. FREE
CONSULTATION
WITH A WATER EXPERT

Order Testing Kits For Accurate Results

- STANDARD KIT (#80101) - Includes mailer kit, lab service, 1 specimen vial & prepaid Fedex overnight shipping. Available in 4 pack (#80401), 6 pack (#80601), and 16 pack (#81601)
- ECONOMY KIT (#70101) - Includes mailer kit, lab service, & 1 specimen vial
- ADDITIONAL VIALS (#80100) - Add as many vials as you need to any of the above standard or economy kits