

## All Gauze Sponges



AVAILABLE IN:

2x2 (8-ply or 12-ply)

3x3 (12-ply)

4x4 (8-ply or 12-ply)

Proof that a product that works well can also be simple and inexpensive. All Gauze Sponges have a loose weave, giving them a textured, cheesecloth-like appearance and a mildly abrasive surface, which means they provide a good grip.

- Simple yet effective
- For use in a wide variety of applications
- Mildly abrasive texture ideal for cleaning instruments and surfaces
- Autoclavable, leak proof dispensers available for 2"x2" and 4"x4" sizes, and come in two colors: blue and white



CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
□ 300420	2"x2"	ALL-GAUZE 12-PLY - NON-STERILE	200 SPONGES/SLEEVE	25 SLVS/CASE	CHINA
□ 300430	2"x2"	ALL-GAUZE 8-PLY - NON-STERILE	200 SPONGES/SLEEVE	25 SLVS/CASE	CHINA
□ 300421	3"x3"	ALL-GAUZE 12-PLY - NON-STERILE	200 SPONGES/SLEEVE	20 SLVS/CASE	CHINA
□ 300422	4"x4"	ALL-GAUZE 8-PLY - NON-STERILE	200 SPONGES/SLEEVE	20 SLVS/CASE	CHINA
□ 300423	4"x4"	ALL-GAUZE 12-PLY - NON-STERILE	200 SPONGES/SLEEVE	10 SLVS/CASE	CHINA
□ 300101	2"x2" SPONGE DISPENSER	BLUE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 300103	2"x2" SPONGE DISPENSER	WHITE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 300111	4"x4" SPONGE DISPENSER	BLUE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 300113	4"x4" SPONGE DISPENSER	WHITE - LEAKPROOF	1 DISPENSER	N/A	USA

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.



RichmondDentalAndMedical.net

855-841-0743

