

## Cotton-Filled Sponges



AVAILABLE IN:

2x2 (8-Ply)  
(Sterile, Non-sterile)

3x3 (4-Ply)  
(Non-sterile)

4x4 (4-Ply)  
(Non-sterile)

This variety of sponge looks similar to the all-gauze variety with its cheesecloth-like outer layer, but contains an inner layer of 100% cotton. This provides extra absorbency, as well as additional patient comfort. The cotton-filled sponge is dentists' traditional, reliable choice.

- All gauze with an inner layer of 100% cotton
- Cotton center provides greater absorbency
- Autoclavable, leak proof dispensers available for 2x2 and 4x4 sizes, and come in two colors: blue and white
- 2"x2" size available in pre-sterilized SteriPocket®. These sponges are pre-packaged sterile sponges and are ideal to distribute to patients after extractions.



CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
□ 300405	2"x2"	COTTON FILLED 8-PLY - <b>STERILE</b>	200 SPONGES/SLEEVE	25 SLVS	CHINA
□ 300410	2"x2"	COTTON FILLED 8-PLY - NON-STERILE	200 SPONGES/SLEEVE	25 SLVS	CHINA
□ 300311	3"x3"	COTTON FILLED 4-PLY - NON-STERILE	200 SPONGES/SLEEVE	20 SLVS	CHINA
□ 300312	4"x4"	COTTON FILLED 4-PLY - NON-STERILE	200 SPONGES/SLEEVE	10 SLVS	CHINA
□ 300101	2"x2" SPONGE DISPENSER	BLUE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 300103	2"x2" SPONGE DISPENSER	WHITE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 300111	4"x4" SPONGE DISPENSER	BLUE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 300113	4"x4" SPONGE DISPENSER	WHITE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 300408	2"x2", COTTON FILLED SPONGES	2 SPONGES/PK - <b>STERILE</b>	300 PKs/BX	8 BOXES	CHINA

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.