

Earloop Face Mask

LEVEL I

AVAILABLE IN:

Anti-fluid
Anti-fluid, Anti-fog
Anti-fluid, Anti-fog
with shield

COLORS:

Pink
Blue
Lavender

They are as integral to a dental procedure as the most expensive equipment in the room. After all, you can't work without protection from airborne germs and fluids. Fortunately, our earloop face masks give you just such protection without sacrificing breathability. For added protection, we offer a face mask with an eye protection shield.

- 100% latex-free, including the earloops
- Face mask pleats cover nose, mouth and chin for full protection
- 3-ply spunbound fabric provides comfort, filtration and excellent breathability
- Particle filtration efficiency of $\geq 95\%$ at 1.0 micron
- Soft aluminum nosepiece maintains shape when formed to face and is fully enclosed for comfort

CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
□ 400501	ANTI-FLUID	PINK	50 MASKS/BOX	6 BOXES	USA
□ 400502	ANTI-FLUID	BLUE	50 MASKS/BOX	6 BOXES	USA
□ 400509	ANTI-FLUID	LAVENDER	50 MASKS/BOX	6 BOXES	USA
□ 400503	ANTI-FLUID, ANTI-FOG	PINK	50 MASKS/BOX	6 BOXES	USA
□ 400504	ANTI-FLUID, ANTI-FOG	BLUE	50 MASKS/BOX	6 BOXES	USA
□ 400510	ANTI-FLUID, ANTI-FOG	LAVENDER	50 MASKS/BOX	6 BOXES	USA
□ 400505	ANTI-FLUID, ANTI-FOG W/SHIELD	PINK	25 MASKS/BOX	4 BOXES	USA
□ 400506	ANTI-FLUID, ANTI-FOG W/SHIELD	BLUE	25 MASKS/BOX	4 BOXES	USA
□ 400511	ANTI-FLUID, ANTI-FOG W/SHIELD	LAVENDER	25 MASKS/BOX	4 BOXES	USA

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.