

Wrapped Rolls



AVAILABLE IN:

Sterile

Non-sterile

DIAMETER:

Medium

LENGTH:

1-1/2"

Our Wrapped Rolls are a great, reliable alternative among dental rolls. Their soft, nonwoven absorbent casing contains no starch which means they will not stick to mucosa.

- 100% purified medical-grade cotton interior wrapped in nonwoven fabric
- Sealed with adhesive approved by the FDA
- Starch-free; will not stick to mucosa
- Infection control roll dispensers are available; they are fully autoclavable and will keep rolls sterile until use; available in blue and white.



CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
200404	1-1/2"	MED. DIA. (3/8") - NON-STERILE	2000/BOX	12 BOXES	USA
200406	1-1/2"	MED. DIA. (3/8") - STERILE	2000/BOX	12 BOXES	USA
200224	INFECTION CONTROL ROLL DISPENSER	BLUE - AUTOCLAVABLE	1 DISP + 200 ROLLS*	6 DISPENSERS	USA
200424	INFECTION CONTROL ROLL DISPENSER	WHITE - AUTOCLAVABLE	1 DISP + 200 ROLLS*	6 DISPENSERS	USA

*Dispensers are shipped with Braided Rolls only.

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.