

FLEXAPRENE* GREEN

Powder-Free Exam Gloves



NEW TECHNOLOGY. NEW LEVEL OF COMFORT.

Wish you could remove the threat of latex allergies but keep the comfortable feel of latex? You just found that glove.

The new FLEXAPRENE* GREEN Exam Glove is an advanced chloroprene technology that gives you the feel and tactile sensitivity of latex without the risk of Type I allergic reactions to natural rubber latex proteins – and with the protection you'd expect from a synthetic.

FEEL THE DIFFERENCE

- Tactile sensitivity comparable to latex
- Soft to the touch, yet powerful in grip
- Polymer coated for improved donning
- Not made with natural rubber latex
- Powder-free chloroprene technology



From the manufacturer of the
#1 branded nitrile glove
in the dental market.¹

NEW!

FLEXAPRENE* GREEN Powder-Free Exam Gloves

A NEW Chloroprene Technology

**SO COMFORTABLE, WE HAD TO MAKE THEM GREEN -
SO YOU'D REMEMBER TO TAKE THEM OFF!**

Can exam gloves really be that comfortable? The answer is yes - with new FLEXAPRENE* GREEN Exam Gloves.

They're soft, easy to don, and fit so comfortably that you may just forget to take them off!



FLEXAPRENE* GREEN Exam Gloves

Item	Size	Packaging	
44792	X-Small	200 eaches/box	10 boxes/case
44793	Small	200 eaches/box	10 boxes/case
44794	Medium	200 eaches/box	10 boxes/case
44795	Large	200 eaches/box	10 boxes/case
44796	X-Large	180 eaches/box	10 boxes/case

Product Specifications

Gauge Thickness Measurements	MM	MIL
Middle Finger:	.09	3.5
Palm:	.07	2.8
Cuff:	.06	2.4
Average Length:	240mm	(9.5")

Physical Properties

Before Aging	
Tensile Strength:	16 MPa
Ultimate Elongation:	550%
After Aging	
Tensile Strength:	17 MPa
Ultimate Elongation:	450%

Quality Standards

Exceeds current ASTM Standards for critical defects, tensile strength, and puncture resistance.

Passes primary skin irritation testing in accordance with ISO 10993-10(2010).

Manufactured in accordance with Quality System ISO 13485.

AQL - ASTM D5151

Critical Defects Level ≤ 1.5 AQL[†]

[†] Acceptable Quality Level



KNOWLEDGE NETWORK* Accredited Education
In-Service Training and Technical Support
Credentialed Sales Representatives
Tools & Best Practices
Clinical Research
Commitment to Excellence

Call 1-844-HALYARD (1-844-425-9273)
in the United States and Canada.

For more information, please visit
www.halyarddental.com
or call your Halyard Sales Representative.

www.halyardhealth.com



[†]Source: Strategic Data Marketing, 2015 Dental Distributor Data for LAVENDER* branded gloves.

*Registered Trademark or Trademark of Halyard Health, Inc. or its affiliates. © 2016 HYH.

All rights reserved.

IP249 C161419