

BLACK-FIRE*

Powder-Free Nitrile Exam Gloves



THE GLOVE THAT PULLS DOUBLE DUTY

As a first responder, you do more than one job. Now your glove can, too. The new BLACK-FIRE* Nitrile Exam Glove protects you when working with patients, when you need to direct traffic — it can even alert you when there's a dangerous rip or tear.


Reversible

BLACK-FIRE* has both a black side and a high-visibility orange side you can use when directing traffic or working in low-visibility situations.



Breach Detection

Proprietary QUICK CHECK* technology lets you quickly spot breaches in protection caused by cuts and tears.



First responders have to be prepared for any situation. Innovative BLACK-FIRE* gloves from Halyard work right along with you, providing dual-purpose protection you need in a glove.



NEW!

BLACK-FIRE*

Powder-Free Nitrile Exam Gloves

BLACK-FIRE* NITRILE EXAM GLOVES ARE DESIGNED FOR SUPERIOR PROTECTION AND SAFETY.

- NFPA 1999-2013 Certified
- Patented nitrile formulation increases comfort and strength
- Not made with natural rubber latex

Two Gloves In One



BLACK-FIRE* dual-purpose gloves can be reversed and worn with the bright orange side out in situations that call for high-visibility.

QUICK CHECK* technology



BLACK-FIRE* gloves with proprietary QUICK CHECK* technology offer the unique ability to detect rips and tears by showing the high-visibility orange inner layer through the breach.

Caution: These exam gloves do not protect against heat or fire. Do not use in applications involving exposure to fire, flames or other sources of heat.



KNOWLEDGE NETWORK* Accredited Education
In-Service Training and Technical Support
Credentialed Sales Representatives
Tools & Best Practices
Clinical Research
Commitment to Excellence

Call 1-844-HALYARD (1-844-425-9273)
in the United States and Canada.

www.halyardhealth.com



BLACK-FIRE* Nitrile Exam Gloves

Item	Size	Packaging	
44755	X-Small	150 eaches/box	10 boxes/case
44756	Small	150 eaches/box	10 boxes/case
44757	Medium	150 eaches/box	10 boxes/case
44758	Large	150 eaches/box	10 boxes/case
44759	X-Large	140 eaches/box	10 boxes/case

Product Specifications

Gauge Thickness Measurements	MM	MIL
Middle Finger:	.14	5.5
Palm:	.10	3.9
Cuff:	.09	3.5
Average Length:	242mm	(9.5")

Physical Properties

Before Aging	
Tensile Strength:	42 MPa
Ultimate Elongation:	650%
After Aging	
Tensile Strength:	38 MPa
Ultimate Elongation:	550%

Quality Standards

Exceeds current ASTM Standards for critical defects, tensile strength, and puncture resistance.

Passes primary skin irritation testing in accordance with ISO 10993-10(2010).

Manufactured in accordance with Quality System ISO 9001.

AQL - ASTM D5151

Critical Defects Level <1.0 AQL[†]

[†] Acceptable Quality Level



Know a first responder deserving recognition? Honor them through **Fist Bump for Heroes**, a program designed to celebrate those on the front lines.

To learn more about BLACK-FIRE* gloves and Fist Bump for Heroes, [visit Black-Fire.com](http://Black-Fire.com).

Patent Pending

*Registered Trademark or Trademark of Halyard Health, Inc. or its affiliates. © 2016 HYH. All rights reserved.

IP249 C161418