




BLACK-FIRE* Nitrile Exam Gloves

THE GLOVE THAT PULLS DOUBLE DUTY.



BLACK SIDE
BREACH
DETECTION



ORANGE SIDE
REVERSIBLE
HIGH-VISIBILITY

AS A FIRST RESPONDER,
YOU DO MORE THAN ONE JOB.

NOW YOUR GLOVES CAN TOO.

- Proprietary QUICK CHECK* technology reveals an orange inner layer when there's a rip or tear
- Orange side can be used for directing traffic or working in low-visibility situations

To learn more about BLACK-FIRE* gloves, visit Black-Fire.com  

Caution: These exam gloves do not protect against heat or fire.
Do not use in applications involving exposure to fire, flames or other sources of heat.

*Registered Trademark or Trademark of Halyard Health, Inc. or its affiliates.
© 2016 HYH. All rights reserved. IP249 C161472